

Pouches, self seal

03



For use in high vacuum, high temperature steam autoclaves and other sterilisation systems.

Self seal pouches are made from the same materials as standard pouches but with a self sealing adhesive strip.

Simple removal of the backing strip, folding over and firm pressure will ensure a correct seal.

Catalogue No	Code	Width, mm	Length, mm	Pack qty
AUY-160-012B	S1	90	203	4,000
AUY-160-051X	S4	133	254	2,000
AUY-160-071Y	S7	190	330	1,000
AUY-160-120V	S2	90	230	2,000
AUY-160-130S	S3	90	559	1,000
AUY-160-140P	S6	133	356	2,000
AUY-160-160J	S8	305	381	1,000

Bags, polyethylene

33



Catalogue No	Width, mm	Depth, mm	Pack qty
BAJ-300-050J	100	150	100
BAJ-300-070D	125	175	100
BAJ-300-090U	200	250	100
BAJ-300-110R	200	300	50
BAJ-300-130L	250	375	50

Bags, autoclave, storage, resealable, write-on and plain, polyethylene, Minigrip®

33



Ideal for storing stationery, sending samples etc.

Polyethene, resealable bags, clear or with three write-on bars for labelling.

Write-on bags

Catalogue No	Width, mm	Depth, mm	Pack qty
BAJ-340-003N	40	65	100
BAJ-340-005J	60	60	100
BAJ-340-011X	75	85	100
BAJ-340-031Y	90	115	100
BAJ-340-040H	115	115	100
BAJ-340-051C	100	140	100
BAJ-340-060B	140	140	100
BAJ-340-065X	75	190	100
BAJ-340-071T	125	190	100
BAJ-340-091N	190	190	100
BAJ-340-110M	150	230	100
BAJ-340-120J	205	280	100
BAJ-340-130G	230	325	100
BAJ-340-140D	57	76	100
BAJ-340-150A	254	356	100
BAJ-340-160U	279	406	100
BAJ-340-170R	330	457	100

Plain bags

Catalogue No	Width, mm	Depth, mm	Pack qty
BAJ-360-010C	40	65	100
BAJ-360-020W	60	60	100
BAJ-360-040Q	60	75	100
BAJ-360-050N	75	80	100
BAJ-360-060K	90	115	100
BAJ-360-070H	115	115	100
BAJ-360-080E	100	140	100
BAJ-360-090B	140	140	100
BAJ-360-110V	75	190	100
BAJ-360-120S	125	190	100
BAJ-360-130P	190	190	100
BAJ-360-140M	150	230	100
BAJ-360-150J	50	230	100
BAJ-360-160G	205	280	100
BAJ-360-170D	230	325	100
BAJ-360-180A	325	325	100
BAJ-360-190U	255	355	100
BAJ-360-210R	280	405	100
BAJ-360-220X	330	455	100
BAJ-360-230L	380	510	100