



Bags, biohazard

56



a

Non-sterile sturdy disposal bags with biohazard labelling.

Available in standard temperature, 121°C, or high temperature, 134°C versions. All 'tissue box' style dispensers. Opening along width.

Standard temperature (121°C), polyethylene bags

Catalogue No	Alt. No	Width, mm	Depth, mm	Pack qty
BAJ-440-010Y	509	305	660	200
BAJ-440-020F	509L	405	610	200
BAJ-440-026Q	510L	405	610	500
BAJ-440-040W	511	610	810	200
BAJ-440-051R	510	305	660	500

High temperature (134°C), polypropylene bags

Catalogue No	Alt. No	Width, mm	Depth, mm	Pack qty
BAJ-440-030C	511HT	610	810	200
BAJ-440-046K	510HT	305	660	500
BAJ-440-056H	510LHT	405	610	500
BAJ-440-060Q	509HT	305	660	200
BAJ-440-070N	509LHT	405	610	200

Accessories

Catalogue No	Alt. No	Description	Pack qty
BAJ-440-500X	S23C	Non-sterile metal holder for bag size 305mm x 660mm and 405mm x 610mm	1
BAJ-440-510L	S23B	Non-sterile cardboard holder for bag size 610mm x 810mm	10
BAJ-440-520Y	S23E	Non-sterile metal holder for bag size 610mm x 810mm	1

Waste bags

03

Polyethylene sacks, bags and bin liners in a variety of colours and weights.

Catalogue No	Description	Width, mm	Depth, mm	Height, mm	Gauge	Pack qty
WDS-800-010U	Low density PE, black	457	737	990	150/160	200
WDS-800-020R	Low density PE, black	457	737	990	180/200	200
WDS-800-030R	Polyethylene, clear	457	737	990	120/140	200
WDS-800-040L	Sacks, yellow, plain	457	737	990	180/200	200
WDS-800-060F	Sacks, green, plain	457	737	990	140/160	200
WDS-800-080W	50L bin liners, white	330	585	712	50	500