

Bottles, borosilicate glass, laboratory, reagent and media, crimp top

50



Rubber stoppers - 3 leg lyophilisation

Catalogue No	Alt. No	Description	For neck dimensions, mm	Pack qty
BTF-784-130G	224100-202	Grey chlorobutyl/50	13 x 20 (I.D. x O.D.)	1,000

Aluminium crimp seals

Catalogue No	Alt. No	Description	For neck dimensions, mm	Pack qty
BTF-785-050V	24178-01	Open top, natural	20 (O.D.)	1,000
BTF-785-070P	224183-01	Centre disc tear out, natural	20 (O.D.)	1,000
BTF-785-090J	224191	Solid top, natural	20 (O.D.)	1,000
BTF-785-110G	224193-01	Tear off, natural	20 (O.D.)	1,000



Bottles, borosilicate, screw neck, diagnostic

03



a

Ideal for the packaging of diagnostic and lyophilised reagents.

Clear or amber bottles with a wide range of specially designed thin flange stoppers and caps. Conform to EP Type 1 standards. Autoclavable.

Clear

Catalogue No	Alt. No	Capacity, mL	Cap size, mm	Height, mm	O.D., mm	Pack qty
BTF-770-050M	219360	5	20-400	48	21	480
BTF-770-070G	219361	10	20-400	61	25	480
BTF-770-090A	219362	20	20-400	66	32	480
BTF-770-110U	219364	50	20-400	80	42	288

Amber

Catalogue No	Alt. No	Capacity, mL	Cap size, mm	Height, mm	O.D., mm	Pack qty
BTF-772-050A	219370	5	20-400	48	21	480
BTF-772-070R	219371	10	20-400	61	25	480
BTF-772-090L	219372	20	20-400	66	32	480
BTF-772-110Y	219374	50	20-400	80	42	288

Rubber stoppers - lyophilisation

Catalogue No	Alt. No	Stopper material	Neck O.D., mm	Pack qty
BTF-774-050L	224100-203	Grey butyl	20	300

Closures - I-Loc™, white polypropylene with grey stopper

Catalogue No	Alt. No	Cap size, mm	Pack qty
BTF-776-050W	240676-02	20-400	100

Closures - open top cap, white polypropylene

Catalogue No	Alt. No	Cap size, mm	Pack qty
BTF-776-150S	240716-02	20-400	300

Closures - closed top linerless cap, white polypropylene

Catalogue No	Alt. No	Cap size, mm	Pack qty
BTF-776-250X	240706-02	20-400	300

