




Accessories and spare parts for BTF-682, BTF-685, BTF-710, BTF-720 and BTF-730 bottles

01



Autoclavable.

Polypropylene screw caps, with lip seal, autoclavable to 140°C, blue

Catalogue No	Alt. No	Aperture size, mm	ISO thread GL	Pack qty
BTF-690-501W	2923913	-	25	10
BTF-690-502U	2923919	-	32	10
BTF-690-504Q	2923928	-	45	10

Polypropylene pouring rings, autoclavable to 140°C, blue (for screw caps above)

Catalogue No	Alt. No	Aperture size, mm	ISO thread GL	Pack qty
BTF-690-540M	2924219	-	32	10
BTF-690-506M	2924228	-	45	10

PBT screw caps, with PTFE protected seals, autoclavable to 180°C, red

Catalogue No	Alt. No	Aperture size, mm	ISO thread GL	Pack qty
BTF-690-520S	2924013	-	25	10
BTF-690-530P	2924019	-	32	10
BTF-690-508Y	2924028	-	45	10

EFTE pouring rings, autoclavable to 180°C, red (for screw caps above)

Catalogue No	Alt. No	Aperture size, mm	ISO thread GL	Pack qty
BTF-690-550J	2924419	-	32	10
BTF-690-510V	2924428	-	45	10

PBT caps with aperture, autoclavable to 180°C, red

Catalogue No	Alt. No	Aperture size, mm	ISO thread GL	Pack qty
BTF-690-512R	2922708	20	32	10
BTF-690-514N	2922710	34	45	10

Sealing rings for aperture caps. Silicone rubber with a PTFE protective face

Catalogue No	Alt. No	Aperture size, mm	ISO thread GL	Pack qty
BTF-690-516J	2923610	10	32	10
BTF-690-522X	2923614	14	32	10
BTF-690-524K	2923616	16	32	10
BTF-690-526G	2923618	18	32	10
BTF-690-528C	2923826	26	45	10
BTF-690-532L	2923832	32	45	10

Septa seals, silicone (for piercing) for aperture caps

Catalogue No	Alt. No	Aperture size, mm	ISO thread GL	Pack qty
BTF-690-534H	2924608	29	32	100
BTF-690-536D	2924610	42	45	100



Premium cap, GL45 thread

01



Manufactured from thermoplast/duraplast (PTch260), this cap has the highest chemical resistance.

The cap contains no colourants, and has a welded-in silicone seal coated with PTFE on both sides. Matching pouring rings are available and guarantee non-dripping and safe handling. The premium cap has a temperature resistance of -196°C to 260°C and therefore suitable for many sterilisation methods such as autoclaving and dry sterilisation. The premium cap meets all the requirements of both US pharmacopoeia (USP class VI) and European pharmacopoeia (EP) and therefore complies with FDA standards. Autoclavable.

Caps

Catalogue No	Alt. No	ISO thread GL	Pack qty
BTF-690-561E	1088679	45	5

Pouring rings

Catalogue No	Alt. No	ISO thread GL	Pack qty
BTF-690-536A	1088678	45	5