

Bottles, polypropylene, aspirators

02

Azelon

a

Squat form aspirators, for easy placement on shelves or bench top, with moulded in handle for safe and secure handling.

Moulded graduations in litres and moulded caps with drill guides for easy location of probes and feed lines. Stopcock opens and closes in a quarter turn. Autoclavable. Closures and stopcocks can be removed or loosened before autoclaving to avoid accidental implosion.

Catalogue No	Alt. No	Capacity, L	Cap size, mm	Height with cap, mm	Width, mm
ASP-755-010A	BNP05A	5	45	270	220
ASP-755-020U	BNP10A	10	90	360	280
ASP-755-030R	BNP20A	20	90	420	350

Spare parts

Catalogue No	Alt. No	Description
ASP-755-040X	CAP045H	Cap for ASP-755-010A/BTK-610-010B
ASP-755-050L	CGZ85PP	Cap for ASP-755-020U/ASP-755-030R/BTK-610-020V/BTK-610-030S
ASP-755-060Y	SCBNP01	Spare stopcock



Bottles, polypropylene, aspirators, Thermo Scientific Nalgene®, Lowboy™

55

Thermo SCIENTIFIC

a

Excellent chemical resistance with leakproof closure and tap.

Autoclavable. Before autoclaving, closure should be disengaged from the threads.

Catalogue No	Alt. No	Capacity, L	Neck size, mm	Height, mm	O.D., mm	Closure size, mm
ASP-810-030V	2324-0008	8	44	165	368	63
ASP-810-050P	2324-0015	15	44	203	368	63

Spare part

Catalogue No	Alt. No	Description
ASP-704-010S	6422-0010	Tap



Bottles, polycarbonate, aspirators, Thermo Scientific Nalgene®, Clearboy™

55

Thermo SCIENTIFIC

a

Transparent and tough, with polypropylene tap and screw closure. Useful for many applications, including use in cold rooms of sub-zero temperatures.

Graduated in 5L increments. Supplied with tap closure. Leakproof. Autoclavable.

Catalogue No	Alt. No	Capacity, L	Neck size, mm	Dimensions [O.D. x h], mm	Closure size, mm
ASP-700-010T	2317-0020	10	64	249 x 389	83B
ASP-700-020Q	2317-0050	20	64	282 x 513	83B

Spare part

Catalogue No	Alt. No	Description
ASP-704-010S	6422-0010	Tap

