

Bottles jerricans



Bottles, polycarbonate, aspirators, rectangular, Thermo Scientific Nalgene®

55

Thermo
SCIENTIFIC

Useful for large volume media and culture preparation especially where light is required, e.g. algae.

Ideal for storage or dispensing reagents or distilled water. Practical for use in cold rooms even at sub-zero temperatures. Graduated in 5L increments. White polypropylene screw closure. Leakproof.

Catalogue No	Alt. No	Capacity, L	Neck size, mm	Height, mm	O.D., mm	Closure size, mm
ASP-850-030Q	2322-0020	9	86.4	362	216	100
ASP-850-050K	2322-0050	20	86.4	381	218	100

Spare part

Catalogue No	Alt. No	Description
ASP-704-010S	6422-0010	Tap



Bottles, polyethylene, jerricans, narrow neck

02

Fisherbrand®

Integrally moulded handle and tamper evident screw cap.

Catalogue No	Capacity, L	Height, mm	Width, mm	Neck O.D., mm
FB59631	1	165	125	38
FB59632	2.5	205	150	45
FB59633	5	282	189	38
FB59634	10	314	230	61
FB59635	25	470	295	61



Bottles, polyethylene, jerricans, narrow neck

03

Integrally moulded handle and tamper evident screw cap. All sizes except 1L are UN approved.

Catalogue No	Capacity, L	Width, mm	Depth, mm	Height, mm	Neck I.D., mm
BTK-525-010S	1	125	70	165	38
BTK-525-020P	2.5	150	110	205	45
BTK-525-050G	5	189	133	282	38
BTK-525-075N	10	231	197	314	61



Bottles, fluorinated HDPE, translucent, jerricans, wide neck, Thermo Scientific Nalgene®

55

Thermo
SCIENTIFIC

Durable, specially designed for hard use.

Surface fluorinated inside and outside, providing barrier to solvent absorption and penetration. Long integral spout allows easy pouring. With fluorinated polypropylene screw closure and built-in handle. Suitable for acids, alkalis and aggressive organic solvents. Leakproof.

Catalogue No	Alt. No	Capacity, L	Height, mm	Neck I.D., mm	Closure
BTK-710-070V	2242-0025	10	368	37	53B
BTK-710-110M	2242-0050	20	445	40	53B