



Bottles, HDPE, translucent, narrow neck, Thermo Scientific Nalgene®

Thermo
SCIENTIFIC

55

Clear and amber bottles with excellent chemical resistance ideal for storing, shipping and packaging liquids. Leakproof. Good freezer storage to -100°C.

Clear

Catalogue No	Alt. No	Capacity, mL	Neck size, mm	Height, mm	O.D., mm	Closure size, mm	Pack qty
BTK-225-005M	2002-9125	4	9	41	16	13	12
BTK-225-010T	2002-9025	8	13	45	23	20	12
BTK-225-020Q	2002-9050	15	13	56	23	20	12
BTK-225-030N	2002-0001	30	13	61	36	20	12
BTK-225-050H	2002-0002	60	13	84	38	20	12
BTK-225-070B	2002-0004	125	18	102	51	24	12
BTK-225-080V	2002-0006	175	18	125	54	24	12
BTK-225-090S	2002-0008	250	18	135	61	24	12
BTK-225-110P	2002-0016	500	23	170	74	28	12
BTK-225-130J	2002-0032	1,000	28	216	91	38	6

Amber

Catalogue No	Alt. No	Capacity, mL	Neck size, mm	Height, mm	O.D., mm	Closure size, mm	Pack qty
BTK-240-001N	2004-9125	4	9	41	16	13	12
BTK-240-003J	2004-9025	8	13	45	23	20	12
BTK-240-005F	2004-9050	15	13	56	23	20	12
BTK-240-008W	2004-0001	30	13	61	36	20	12
BTK-240-010M	2004-0002	60	13	84	38	20	12
BTK-240-020J	2004-0004	125	18	102	51	24	12
BTK-240-030G	2004-0008	250	18	135	61	24	12
BTK-240-050A	2004-0016	500	23	170	74	28	12
BTK-240-070R	2004-0032	1,000	28	216	91	38	6



Bottles, HDPE, translucent, narrow neck

41

Fisherbrand

General purpose bottles with polypropylene screw closure.

Catalogue No	Alt. No	Capacity, mL	Neck size, mm	Height, mm	O.D., mm	Pack qty
FB59611	P109/C1	30	18/R4	68	30	72
FB59613	P109/C3	125	22/R6	102	48	72
FB59614	P109/C5	250	22/R6	130	59	72
FB59615	P109/C6	500	24/R4	170	73	48
FB59616	P109/C7	1,000	28/R6	248	83	24



Bottles, HDPE, translucent, medium weight, narrow neck

02

Azlon

Ideal for general liquid storage and sampling.

Supplied capped and in shrink-wrap packs. All sizes feature S.P.I. recycling codes on their base. Polypropylene screw closure. Leakproof.

Catalogue No	Alt. No	Capacity, mL	Neck size, mm	Height, mm	O.D., mm	Pack qty
BTK-285-030R	BNH0030PN	30	16	67	32	10
BTK-285-050L	BNH0060PN	60	16	78	40	10
BTK-285-070F	BNH0125PN	125	21	103	50	10
BTK-285-090W	BNH0150N	150	21	116	50	10
BTK-285-110T	BNH0250P	250	21	140	60	10
BTK-285-130N	BNH0500P	500	21	160	75	10
BTK-285-150H	BNH1000P	1,000	21	201	90	5
BTK-285-170B	BNH2000P	2,000	34	237	120	5