

Bottles, LDPE, translucent, wide neck

41

Fisherbrand®

Clear bottles with blue polypropylene screw closure. Leakproof.

Catalogue No	Capacity, mL	Neck size, mm	Height, mm	O.D., mm	Pack qty
FB59636	30	33/R3	45	36	12
FB59637	60	33/R3	80	36	12
FB59638	125	38/R3	93	50	12
FB59639	250	58/R3	114	61	12
FB59640	500	58/R3	140	76	12
FB59641	1,000	70/R3	200	90	6



Bottles, LDPE, translucent, wide neck, Thermo Scientific Nalgene®

55

**Thermo
SCIENTIFIC**

Flexible with excellent impact resistance, ideal for shipping liquids. LDPE and polypropylene screw closure. Leakproof.

Catalogue No	Alt. No	Capacity, mL	Neck size, mm	Height, mm	O.D., mm	Closure size, mm	Pack qty
BTK-345-030X	2103-0001	30	23	64	36	28	12
BTK-345-050Y	2103-0002	60	20	86	38	28	12
BTK-345-070C	2103-0004	125	31	99	51	38	12
BTK-345-090T	2103-0008	250	33	130	61	43	12
BTK-345-110Q	2103-0016	500	43	168	74	53	12
BTK-345-130K	2103-0032	1,000	53	201	91	63	6



For sampling scoops see page 1,252 - 1,254

Bottles, LDPE, translucent, wide neck

02

Kartell

For solid and liquid samples.

Withstand considerable pressure before rupture. Used by horse racing laboratories for transporting sensitive liquid and solid samples from course to laboratory. Secured by tagging points on bottle and cap. Confidential notes can be placed in the space between insert and sealing plug and cap. Graduated.

Catalogue No	Alt. No	Capacity, mL	Neck size, mm	Height, mm	O.D., mm	Pack qty
BTK-365-010G	408	50	24	75	40	10
BTK-365-020D	409	100	24	90	50	10
BTK-365-030A	410	250	34	130	60	10
BTK-365-040U	411	500	34	160	75	10
BTK-365-051P	412	1,000	44	200	95	5
BTK-365-061M	413	2,000	44	250	115	1



Bottles, LDPE, translucent, wide neck

02

Kartell

Graduated general purpose bottles. Long pinch seal ensures leakproof closure.

Conforms to DIN 13316 and 168.

Catalogue No	Alt. No	Capacity, mL	Neck size, mm	Height, mm	O.D., mm	Pack qty
BTK-385-010P	1608	50	24	88	38	10
BTK-385-020M	1610	100	24	103	48	10
BTK-385-030J	1612	250	38	140	60	10
BTK-385-040G	1614	500	38	168	75	10
BTK-385-051B	1616	1,000	55	206	95	5
BTK-385-061V	1618	2,000	55	249	120	1

