

Bottles HDPE



Bottles, HDPE, translucent, wide neck, Thermo Scientific Nalgene®

55

Thermo
SCIENTIFIC

Clear and amber durable, general purpose bottles.
Suitable for shipping liquids. Polypropylene screw closure.

Clear

Catalogue No	Alt. No	Capacity, mL	Neck I.D., mm	Height, mm	O.D., mm	Closure size, mm	Pack qty
BTK-478-030N	2104-0001	30	23	63	36	28	12
BTK-478-050H	2104-0002	60	20	86	38	28	12
BTK-478-070B	2104-0004	125	31	99	51	38	12
BTK-478-090S	2104-0008	250	33	131	61	43	12
BTK-478-110P	2104-0016	500	43	168	74	53	12
BTK-478-130J	2104-0032	1,000	53	199	91	63	6

Amber

Catalogue No	Alt. No	Capacity, mL	Neck I.D., mm	Height, mm	O.D., mm	Closure size, mm	Pack qty
BTK-480-010X	2106-0002	60	20	86	39	28	12
BTK-480-020L	2106-0004	125	31	99	50	38	12
BTK-480-030Y	2106-0008	250	33	131	61	43	12
BTK-480-040F	2106-0016	500	43	168	73	53	12
BTK-480-070T	2106-0032	1,000	53	199	91	63	6



Bottles, HDPE, translucent, wide neck

41

Fisherbrand

Clear general purpose bottles.
Polypropylene screw closure.

Catalogue No	Capacity, mL	Neck size, mm	Height, mm	O.D., mm	Pack qty
FB59602	125	38/R3	93	50	12
FB59603	250	58/R3	114	61	12
FB59604	500	58/R3	140	76	12
FB59605	1,000	70/R3	200	90	6



Bottles, fluorinated HDPE, translucent, wide neck, Thermo Scientific Nalgene®

55

Thermo
SCIENTIFIC

Fluorinated surface provides improved barrier properties and reduces solvent absorption and penetration. Fluorination enhances long term container performance and prevents material loss due to permeation. Split and puncture resistant. Can be used with most acids, alkalis and aggressive organic solvents. Leakproof.

Catalogue No	Alt. No	Capacity, mL	Neck size, mm	Height, mm	O.D., mm	Closure size, mm	Pack qty
BTR-660-030F	2197-0004	125	28	100	50	38	12
BTR-660-050W	2197-0008	250	33	131	61	43	12
BTR-660-070Q	2197-0016	500	43	168	73	53	12
BTR-660-090K	2197-0032	1,000	53	199	92	63	6



Bottles, HDPE, translucent, square, wide neck, Thermo Scientific Nalgene®

55

Thermo
SCIENTIFIC

Square shape saves space.
Wide mouth easy to fill. Leakproof.

Catalogue No	Alt. No	Capacity, mL	Neck size, mm	Height, mm	O.D., mm	Closure size, mm	Pack qty
BTR-310-050M	2114-0002	60	23	81	36	28	12
BTR-310-070G	2114-0006	175	28	102	51	38	12
BTR-310-090A	2114-0008	250	33	114	61	43	12
BTR-310-110U	2114-0016	500	43	140	74	53	12
BTR-310-130X	2114-0032	1,000	53	175	91	63	6