

## Bottles, HDPE, single use, wide neck, round and square, with assembled cap



- Single use, HDPE, round and square bottles from 50mL to 2,500mL
- Available with a choice of three different caps:
  - Single closure, red, deep thread screw cap
  - Double closure, red, polyethylene seal and screw cap
  - Leak proof screw cap, red, with a tresylene seal (100mL to 1,000mL sizes only)
- Bottles are screw cap assembled and set in cases
- Available non-sterile and sterile by irradiation with peel away sterility certificate
- 2,500mL bottles have white caps
- For research use only, unless otherwise indicated

## Technical Specification - Specific

## Round bottle

Capacity, mL	Dimensions, mm	Neck I.D., mm
50	70 x 39	27
100	80 x 52	37
150	108 x 52	37
250	125 x 61	37
500	136 x 81	58
1,000	174 x 100	58
2,000	220 x 122	57

## Round bottle, single closure screw cap, natural

Catalogue No	Capacity, mL	Pack qty
FB55732	100	335
FB55733	150	250
FB55734	250	145
FB55735	500	140
FB55736	1,000	68
FB55737	2,000	25

## Round bottle, double closure: polyethylene seal and screw cap

Catalogue No	Capacity, mL	Pack qty
FB55741	50	450
FB55742	100	335
FB55743	150	250
FB55744	250	145
FB55745	500	140
FB55746	1,000	68
FB55747	2,000	25

## Round bottle, leakproof screw cap with tresylene seal

Catalogue No	Capacity, mL	Pack qty
FB55752	100	335
FB55753	150	250
FB55754	250	145
FB55755	500	140
FB55756	1,000	68

Round bottle - sterile by irradiation with peel away sterilisation certificate  
Single closure screw cap

Catalogue No	Capacity, mL	Pack qty
FB55762	100	335
FB55763	150	250
FB55764	250	145
FB55765	500	140
FB55766	1,000	68
FB55767	2,000	25



## Square bottle

Capacity, mL	Dimensions [w x d x h], mm	Neck I.D., mm
150	50 x 50 x 95	37
250	58 x 58 x 110	37
500	76 x 76 x 122	58
1,000	95 x 76 x 185	58
2,000	109 x 96 x 230	58
2,500	110 x 110 x 265	85



## Round bottle, double closure: polyethylene seal and screw cap

Catalogue No	Capacity, mL	Cap size	Pack qty
FB55771	50	wide	450
FB55772	100	narrow	335
FB55773	150	narrow	250
FB55774	250	narrow	145
FB55775	500	narrow	140
FB55776	1,000	narrow	68
FB55777	2,000	wide	25

## Round bottle - sterile by irradiation with peel away sterilisation certificate

Catalogue No	Capacity, mL	Pack qty
FB55782	100	335
FB55783	150	250
FB55784	250	145
FB55785	500	140
FB55786	1,000	68

## Square, wide mouthed bottle, natural, with single closure leakproof screw cap with tresylene seal.

Catalogue No	Capacity, mL	Pack qty
FB55791	150	300
FB55792	250	210
FB55793	500	175
FB55794	1,000	90
FB55795	2,000	50

entry continued