

Bottles polypropylene



Bottles, polypropylene, translucent, wide neck,
Thermo Scientific Nalgene®

55

Thermo
SCIENTIFIC



For chemicals, specimens and general use.

Leakproof. Autoclavable with cap loosened off threads.

Catalogue No	Alt. No	Capacity, mL	Neck size, mm	Height, mm	O.D., mm	Closure size, mm	Pack qty
BTR-635-030T	2105-0001	30	23	64	36	28	12
BTR-635-050N	2105-0002	60	20	86	38	28	12
BTR-635-070H	2105-0004	125	31	99	51	38	12
BTR-635-090B	2105-0008	250	33	130	61	43	12
BTR-635-110V	2105-0016	500	43	168	74	53	12
BTR-635-130P	2105-0032	1,000	53	201	91	63	6



Bottles, polypropylene, translucent,
heavy weight, wide neck

41

Fisherbrand



Ideal for general laboratory applications.

Polypropylene screw closure. Leakproof. Autoclavable with cap loosened. Excellent chemical resistance.

Clear

Catalogue No	Capacity, mL	Cap size, mm	Height, mm	O.D., mm	Pack qty
FB59642	30	33/R3	45	36	10
FB59643	60	33/R3	80	36	10
FB59644	125	38/R3	93	50	10
FB59645	150	58/R3	68	60	10
FB59646	250	58/R3	114	61	10
FB59647	500	58/R3	140	76	10
FB59648	1,000	70/R3	200	90	5



Bottles, polypropylene, translucent, wide neck

02

Azlon



Wide mouth bottles ideal for liquids and solids.

Manufactured from high clarity polypropylene with good chemical resistance. Clear graduations showing mL and fl oz.

Large 'write-on' space for safe easy labelling, S.P.I. recycling codes moulded to base. Leakproof. Autoclavable to 121°C

Catalogue No	Alt. No	Capacity, mL	Height, mm	O.D., mm	Neck I.D., mm	Pack qty
BTK-396-010C	BWP0060PGRAD	60	82	40	26	10
BTK-396-020C	BLP0125PGRAD	125	100	50	32	10
BTK-396-030T	BWP025PGRAD	250	141	60	34	10
BTK-396-040Q	BWP0060PGAD	500	161	75	44	10
BTK-396-050N	BWP1000GRAD	1,000	207	90	53	5