

## Bottles jars



Jars, polypropylene, wide neck, Thermo Scientific Nalgene®

**Thermo**  
SCIENTIFIC

55

**a**

| Catalogue No | Alt. No   | Capacity, mL | Dimensions [O.D. x h], mm | Neck I.D., mm | Pack qty |
|--------------|-----------|--------------|---------------------------|---------------|----------|
| JAR-660-010J | 2118-0001 | 30           | 36 x 48 (with closure)    | 36            | 12       |
| JAR-660-030D | 2118-0002 | 60           | 48 x 48 (with closure)    | 43            | 12       |
| JAR-660-050U | 2118-0004 | 125          | 76 x 74 (with closure)    | 64            | 12       |
| JAR-660-070X | 2118-0008 | 250          | 66 x 119 (with closure)   | 64            | 6        |
| JAR-660-090Y | 2118-0016 | 500          | 117 x 86 (with closure)   | 112           | 6        |
| JAR-660-110F | 2118-0032 | 1,000        | 117 x 145 (with closure)  | 112           | 6        |



Bottles, polystyrene, wide neck

**Sterilin**

02

White urea caps with steran faced woodpulp wad.

| Catalogue No | Capacity, mL | Neck size, mm | Height, mm | O.D., mm | Pack qty |
|--------------|--------------|---------------|------------|----------|----------|
| BTK-352-030C | 30           | 26            | 48         | 35       | 20       |
| BTK-352-050T | 60           | 32            | 59         | 43       | 20       |
| BTK-352-070N | 125          | 36            | 68         | 54       | 20       |
| BTK-352-090H | 250          | 54            | 79         | 69       | 10       |
| BTK-352-110E | 350          | 64            | 87         | 80       | 10       |



Jars, PMP, wide neck, Thermo Scientific Nalgene®

**Thermo**  
SCIENTIFIC

55

**a**

Straight sided jars offering good chemical resistance with polypropylene screw closure. Leakproof. Autoclavable, completely disengage threads or remove closure before autoclaving.

| Catalogue No | Alt. No   | Capacity, mL | Height, mm | O.D., mm | Pack qty |
|--------------|-----------|--------------|------------|----------|----------|
| JAR-650-006H | 2117-0060 | 60           | 48         | 48       | 4        |

  

| Catalogue No | Alt. No   | Capacity, mL | Height, mm | O.D., mm | Pack qty |
|--------------|-----------|--------------|------------|----------|----------|
| JAR-650-010Q | 2117-0125 | 125          | 74         | 69       | 4        |
| JAR-650-020N | 2117-0250 | 250          | 119        | 66       | 4        |

  

| Catalogue No | Alt. No   | Capacity, mL | Height, mm | O.D., mm | Pack qty |
|--------------|-----------|--------------|------------|----------|----------|
| JAR-650-030K | 2117-0500 | 500          | 86         | 117      | 4        |
| JAR-650-040H | 2117-1000 | 1,000        | 145        | 117      | 4        |



Bottles, PVC, square, wide neck

02

**Fisherbrand**

Excellent for both organic and inorganic acids, oxidising agents, aldehydes, amines, aliphatic hydrocarbons and mineral oils. Especially useful for powders, pastes and solids.

Generous wide neck bottle with polypropylene cap. Lightweight. Leakproof.

| Catalogue No | Alt. No | Capacity, mL | Neck size, mm | Height, mm | O.D., mm  |
|--------------|---------|--------------|---------------|------------|-----------|
| BTK-340-070J | 1258C   | 100          | 31            | 74         | 46 x 46   |
| BTK-340-090D | 1259C   | 300          | 41            | 98         | 66 x 66   |
| BTK-340-110A | 1260C   | 500          | 54            | 112        | 80 x 80   |
| BTK-340-130R | 1261C   | 1,000        | 69            | 145        | 95 x 95   |
| BTK-340-150L | 1262C   | 2,000        | 69            | 190        | 120 x 120 |