

## Bottle, polyethylene, wash, trigger operated

02

With adjustable nozzle which produces a powerful intermittent or continuous jet or fine spray, head is mounted on a polyethylene bottle shaped to fit the hand.  
Delivers 1.4mL per stroke.

Catalogue No	Capacity, mL	Bottle colour	Height, mm	Width, mm	Depth, mm
WBS-520-J	500	Natural	206	90 (overall)	54



## Bottles, LDPE, wash, oval, narrow neck

02

Requires only slight pressure to dispense liquid.  
For general purpose use.

Catalogue No	Alt. No	Capacity, mL	Pack qty
BTS-580-030M	L780/01	250	5
BTS-580-050G	L780/02	500	5



## Bottles, LDPE, wash, narrow neck

02

Ergonomic wash bottle with integral precision dispensing tube and removable jet tip.  
Suitable for small hands.

Catalogue No	Alt. No	Capacity, mL	Bottle colour	Pack qty
BTS-680-030E	BGF130P	250	Natural	5
BTS-680-210C	BGF136P	500	Natural	5

*Azlon*



## Bottles, LDPE, wash, safety labelled, wide mouthed, colour coded

41

Solvent dispensing

- Leakproof dispensing of your most common solvents
- Easy to read LDPE solvent-resistant printed bottles
- Colour coded for easy identification

Catalogue No	Capacity, mL	Wording	Closure colour	Pack qty
FB56781	500	Distilled Water	Blue	6
FB56783	500	Water	Blue	6
FB56785	1,000	Acetone	Red	6
FB56787	500	Acetone	Red	6
FB56789	500	Methanol	Green	6
FB56791	500	Isopropanol	Yellow	6
FB56793	500	Toluene	Red	6
FB56795	500	Deionised water	Blue	6
FB56797	500	Ethanol	Green	6

**Fisherbrand®**

