

2^{yr}

CFA-152-010N - Swing-bucket rotor (A-4-62-MTP)

Max. rotor capacity and tube types	Dimensions (l x w x h), mm	Max. speed, rpm	Max. RCF, xg	Bucket required
16 x microtitre plates	128 x 85.6 x 15	4,000	2,750	CFA-152-500T/CFA-152-504L
8 x cell culture plates	128 x 85.6 x 21	4,000	2,750	CFA-152-500T/CFA-152-504L
4 x deep well plates	128 x 85.6 x 45	4,000	2,750	CFA-152-500T/CFA-152-504L

Buckets for CFA-152-010N (A-4-62-MTP)

Catalogue No	Alt. No	Description	Pack qty
CFA-152-500T	5810702003	Microtitre plate bucket	4
CFA-152-504L	5810703000	Microtitre plate bucket	1

CFA-146-010T - Swing-bucket rotor (A-4-44)

Max. rotor capacity and tube types	Max. tube size, dia., mm	Max. tube size, length, mm	Max. speed, rpm	Max. RCF, xg	Bucket required	Adapter required
56 x 1.2 to 5.0mL tubes	11	105	5,000	4,200	CFA-145-500Y	CFA-146-504R
48 x 1.5 to 2.0mL tubes	11	111	5,000	4,200	CFA-145-500Y	CFA-146-502V



1 Buckets, caps and adapters for CFA-146-010T (A-4-44)

Catalogue No	Alt. No	Description	Pack qty
CFA-145-500Y	5804741005	Buckets, 100mL to fit A-4-44	4
CFA-146-500C	5804712005	Aerosol-tight caps for 100mL buckets	1
CFA-146-502V	5804751000	Adapter for 1.5 to 2mL tubes	1
CFA-146-504R	5804750004	Adapter for 1.2 to 5mL tubes	1
CFA-146-506N	5804752007	Adapter for 2.6 to 7mL tubes	1
CFA-146-508J	5804753003	Adapter for 3 to 15mL tubes	1
CFA-146-510W	5804717007	Adapter for 15mL Falcon® tubes	1
CFA-146-512S	5804755006	Adapter for 15mL Falcon® tubes	1
CFA-146-514X	5804754000	Adapter for 7 to 17mL tubes	1
CFA-146-516K	5804756002	Adapter for 7 to 18mL tubes	1
CFA-146-518G	5804757009	Adapter for 12 to 30mL tubes	1
CFA-146-520T	5804759001	Adapter for 30 to 50mL tubes	1
CFA-146-522P	5804760000	Adapter for 50 to 75mL tubes	1
CFA-146-524L	5804718003	Adapter for 50mL Falcon® tubes	1
CFA-146-526H	5804758005	Adapter for 50mL Falcon® tubes	1
CFA-146-528D	5804761006	Adapter for 80 to 120mL tubes	1
CFA-146-530Q	5804706005	Bucket, for 2.5mL Falcon® tubes, set of 4	1

CFA-154-010B - Swing-bucket rotor (A-2-DWP)

Max rotor capacity and tube types	Dimensions (l x w x h), mm	Max. speed, rpm	Max. RCF, xg	Bucket required	Adapter required
10 x MTP	128 x 85.6 x 15	3,700	2,250	CFA-154-500H	-
8 x cell culture plates	128 x 85.6 x 21	3,700	2,250	CFA-154-500H	-
4 x deep well plates	128 x 85.6 x 45	3,700	2,250	CFA-154-500H	-
2 x kits	128 x 85.6 x 60	3,700	2,250	CFA-154-500H	-
4 x 384 well PCR* plates	-	-	2,700	CFA-152-527W/CFA-154-506S	CFA-152-500T
4 x 96 well PCR* plates	-	-	2,600	CFA-152-527W/CFA-154-508X	CFA-152-500T
48 x CombiSlides®	-	-	1,000	CFA-152-527W/CFA-154-510E	CFA-152-500T

Buckets, caps and adapters for CFA-154-010B (A-2-DWP)

Catalogue No	Alt. No	Description	Pack qty
CFA-154-500H	5804743008	Deep well plate buckets for A-2-DWP	2
CFA-154-506S	5825713001	Adapter for 384 well PCR* plate	2
CFA-154-508X	5825711009	Adapter for 96 well PCR* plate	2
CFA-154-510E	5825706005	Adapter for CombiSlide® adapter	2
CFA-152-500T	5810702003	Microtitre plate bucket	1
CFA-152-527W	5810741009	Microtitre plate flex bucket	4

CFA-158-010A - Drum rotor (T-60-11)

Max rotor capacity and tube types	Max. tube size, dia., mm	Max. tube size, length, mm	Max. speed, rpm	Max. RCF, xg	Adapter required
120 x 0.4mL tubes	6	120	14,000	16,400	CFA-140-530G
60 x 1.5 to 2mL tubes	11	120	14,000	16,400	CFA-140-528Q



2 Adapters for CFA-158-010A (T-60-11)

Catalogue No	Alt. No	Description	Pack qty
CFA-140-528Q	5804731000	Adapter for 10 x 1.5/2mL test tubes	1
CFA-140-530G	5804732006	Adapter for 20 x 0.4mL test tubes	1