


Centrifuge tubes

Tube, centrifuge

a



Polypropylene - with flat caps

66

CORNING

Catalogue No	Alt. No	Diagram number	Capacity [nominal], mL	Graduations, mL	Dimensions [O.D. x h], mm	RCF, x g [max.]	Pack qty
CFT-420-051C	430828	9	50	0 to 50 x 5	29 x 116 (approx.)	15,500	500
CFT-645-010C	430829	9	50	0 to 50 x 5	35 x 117 (approx.)	15,500	500
TUL-165-090X	4558	9	50	0 to 50 x 5	-	15,500	300

Polypropylene - with plug seal

Catalogue No	Alt. No	Diagram number	Capacity [nominal], mL	Graduations, mL	Dimensions [O.D. x h], mm	RCF, x g [max.]	Pack qty
CFT-643-021J	430290	9	50	0 to 50 x 5	-	15,500	500
CFT-643-011M	430291	9	50	0 to 50 x 5	-	15,500	500

PET - with plug seals

Catalogue No	Alt. No	Diagram number	Capacity [nominal], mL	Graduations, mL	Dimensions [O.D. x h], mm	RCF, x g [max.]	Pack qty
CFT-647-011L	430304	9	50	-	-	3,600	500

Conical. Double threaded caps 15mL and 50mL only.

02

Soda glass - without caps

Catalogue No	Alt. No	Diagram number	Capacity [nominal], mL	Graduations, mL	Dimensions [O.D. x h], mm	RCF, x g [max.]	Pack qty
CFT-390-B	-	2	13	-	17 x 111 (approx.)	1,300	10
CFT-480-A	-	3	15	-	17 x 111 (approx.)	1,300	10
CFT-530-E	-	2	15	0 to 10 x 0.1	17 x 111 (approx.)	1,300	10

Autoclavable. Polymethylpentene can withstand temperatures of up to 170°C. High chemical resistance. Cannot be used with naked flame. Can also be used as conventional test tubes. Crystal clear.

02

Polymethylpentene - without caps

Catalogue No	Alt. No	Diagram number	Capacity [nominal], mL	Graduations, mL	Dimensions [O.D. x h], mm	RCF, x g [max.]	Pack qty
CFT-234-030S	1387	2	10	0.1	16 x 110 (approx.)	-	10
CFT-234-050M	1388	2	15	0.2	18 x 117 (approx.)	-	10

Polypropylene copolymer - with caps

Catalogue No	Alt. No	Diagram number	Capacity [nominal], mL	Graduations, mL	Dimensions [O.D. x h], mm	RCF, x g [max.]	Pack qty
CFS-130-A	3120-0500	7	500	-	73 x 169 (approx.)	4,800	4
CFT-892-S	3120-0250	7	250	-	62 x 128 (approx.)	up to 13,200	4
CFT-920-C	3120-9500	7	500	-	70 x 170 (approx.)	4,800	4
CFT-960-U	3120-1000	7	1,000	-	96 x 184 (approx.)	7,100	4

Polypropylene - without caps

Catalogue No	Alt. No	Diagram number	Capacity [nominal], mL	Graduations, mL	Dimensions [O.D. x h], mm	RCF, x g [max.]	Pack qty
CFT-220-030D	302	3	7	-	12 x 100 (approx.)	-	10
CFT-220-070X	303	3	16	-	18 x 100 (approx.)	-	10
CFT-220-090Y	305	3	31	-	24 x 93 (approx.)	-	10
CFT-220-110F	306	3	55	-	30 x 100 (approx.)	-	10
CFT-220-130W	307	3	75	-	35 x 100 (approx.)	-	10
CFT-220-150Q	308	3	115	-	40 x 120 (approx.)	-	10
CFT-224-030C	2300	2	10	0.1	16 x 110 (approx.)	-	10
CFT-224-050T	2302	2	15	0.2	18 x 117 (approx.)	-	10
CFT-340-N	300	2	10	-	16 x 108 (approx.)	1,300	10