


Tubes, centrifuge

a



Pyrex® glass - without caps

01

PYREX®

Catalogue No	Alt. No	Diagram number	Capacity [nominal], mL	Graduations, mL	Dimensions [O.D. x h], mm	RCF, x g [max.]	Pack qty
CFS-104-U	1660/02M	1	10	-	16 x 90 (approx.)	6,200	1
CFS-120-H	1660/06M	1	25	-	24 x 90 (approx.)	3,000	1
CFT-590-Y	3460/02M	2	15	Yes	17 x 110 (approx.)	3,000	1
CFT-970-N	-	8	10	0 to 0.5 x 0.01	17.5 x 117 (approx.)	1,300	1

The polymethylpentene tubes are autoclavable and can withstand temperatures of up to 170°C. High chemical resistance. Cannot be used with naked flame. Can also be used as conventional test tubes. Crystal clear.

55

Thermo
SCIENTIFIC

Polypropylene copolymer - without caps

Catalogue No	Alt. No	Diagram number	Capacity [nominal], mL	Graduations, mL	Dimensions [O.D. x h], mm	RCF, x g [max.]	Pack qty
CFS-101-040D	3110-0120	1	12	-	16 x 101 (approx.)	50,000	10
CFS-101-050A	3110-0150	1	15	-	16 x 116 (approx.)	50,000	10
CFS-101-090L	3110-0500	1	50	-	29 x 103 (approx.)	50,000	10
CFS-101-100L	3110-1000	1	100	-	32 x 166 (approx.)	50,000	10
CFT-430-M	3103-0015	7	15	Yes	17 x 119 (approx.)	6,000	10
CFT-650-F	3103-0050	7	50	Yes	29 x 134 (approx.)	6,000	10

Polypropylene co-polymer Oak Ridge, tubes - with reusable, linerless non-contaminating 4mm screw caps

Catalogue No	Alt. No	Diagram number	Capacity [nominal], mL	Graduations, mL	Dimensions [O.D. x h], mm	RCF, x g [max.]	Pack qty
CFT-110-050M	3148-0050	9	35	5 to 35 x 1	28.5 x 113 (approx.)	50,000	10

Adapter

Catalogue No	Alt. No	Diagram number	Capacity [nominal], mL	Graduations, mL	Dimensions [O.D. x h], mm	RCF, x g [max.]	Pack qty
CFT-114-510D	DS3147-0050	14	-	-	-	-	6

Polycarbonate, conical - without caps

Catalogue No	Alt. No	Diagram number	Capacity [nominal], mL	Graduations, mL	Dimensions [O.D. x h], mm	RCF, x g [max.]	Pack qty
CFT-440-F	3105-0015	5	15	Yes	17 x 119 (approx.)	6,000	10
CFT-654-E	3105-0050	5	50	Yes	29 x 134 (approx.)	6,000	10

Polycarbonate, round - without caps

Catalogue No	Alt. No	Diagram number	Capacity [nominal], mL	Graduations, mL	Dimensions [O.D. x h], mm	RCF, x g [max.]	Pack qty
CFS-106-Y	3117-0120	1	12	-	16 x 101 (approx.)	50,000	10
CFT-450-V	3117-0150	1	15	-	16 x 116 (approx.)	50,000	10
CFT-631-N	3117-9500	1	50	-	29 x 105 (approx.)	50,000	10
CFT-635-M	3117-0500	1	50	-	29 x 103 (approx.)	50,000	10

Bottle, polycarbonate, polypropylene - with linerless caps

Catalogue No	Alt. No	Diagram number	Capacity [nominal], mL	Graduations, mL	Dimensions [O.D. x h], mm	RCF, x g [max.]	Pack qty
CFT-891-V	3122-0250	7	250	-	62 x 127 (approx.)	27,500	4
CFT-940-L	3122-0500	7	500	-	70 x 168 (approx.)	13,700	4
CFT-980-G	3122-1000	7	1,000	-	96 x 188 (approx.)	7,100	4

Oak Ridge, polycarbonate - with polypropylene linerless caps

Catalogue No	Alt. No	Diagram number	Capacity [nominal], mL	Graduations, mL	Dimensions [O.D. x h], mm	RCF, x g [max.]	Pack qty
CFT-595-010V	3118-0010	4	10	-	16 x 80 (approx.)	50,000	10
CFT-595-030P	3118-0030	4	30	-	26 x 92 (approx.)	50,000	10
CFT-595-040M	3118-0050	4	43	-	29 x 114 (approx.)	50,000	10

Oak Ridge, tubes, polycarbonate - with polypropylene caps

Catalogue No	Alt. No	Diagram number	Capacity [nominal], mL	Graduations, mL	Dimensions [O.D. x h], mm	RCF, x g [max.]	Pack qty
CFT-110-030S	3146-0050	9	35	Yes	28.5 x 104 (approx.)	50,000	10