

Centrifuge tubes



Tubes, centrifuge, 15mL


66

CORNING

S

- Black printed graduations and a large white marking area
- Certified non-pyrogenic and RNase/DNase free
- Sterile by gamma irradiation
- Available in racks or bulk packed

PET = Polyethylene terephthalate
PP = Polypropylene



Catalogue No	Alt. No	Material	Cap style	RCF, x g [max.]	Pack type	Inner pack qty	Pack qty
CFT-415-010Q	430055	PET	Plug seal	3,600	Rack	50	500
CFT-420-070V	430053	PET	Plug seal	3,600	Bag	25	500
CFT-420-075L	430790	PP	CentriStar™ cap	12,000	Rack	50	500
CFT-420-077H	430791	PP	CentriStar™ cap	12,000	Bag	25	500
CFT-420-079D	430766	PP	Plug seal	12,000	Bag	25	500
CFT-420-031Y	430052	PP	Plug seal	12,000	Rack	50	500

Accessories

Catalogue No	Alt. No	Description	Pack qty
TUL-165-520P	431355	Standard rack, 15mL, polystyrene foam	20
TUL-165-510S	4366	Universal rack, 15mL/50mL, polystyrene foam	20



Tubes, centrifuge, 50mL

66

CORNING

S

- Black printed graduations and large white marking spot
- Choice of cap styles
- Original plug seal or flat top rim seal
- Plug seal caps can be colour-coded using inserts
- Available in racks or bulk packed
- Certified non-pyrogenic and RNase/DNase free
- Sterile by gamma irradiation

PET = Polyethylene terephthalate
PP = Polypropylene

Catalogue No	Alt. No	Material	Cap style	RCF, x g [max.]	Pack type	Inner pack qty	Pack qty
CFT-645-010C	430829	PP	CentriStar™	15,500	Bag	25	500
CFT-643-021J	430290	PP	Plug seal	15,500	Rack	25	500
CFT-643-011M	430291	PP	Plug seal	15,500	Bag	25	500
CFT-647-011L	430304	PET	Plug seal	3,600	Rack	25	500
CFT-420-051C	430828	PP	CentriStar™	15,500	Rack	25	500
TUL-165-090X	4558	PP	CentriStar™	15,500	Rack	25	300

Accessories

Catalogue No	Alt. No	Description	Pack qty
TUL-165-500V	4365	Standard rack, 50mL, polystyrene foam	20
TUL-165-510S	4366	Universal rack, 15mL/50mL, polystyrene foam	20