

## Tubes, centrifuge, 50mL, self-standing

66 CORNING

### S

- Black printed graduations with large white marking area
- Certified non-pyrogenic, sterile and RNase/DNase free
- Sterile by gamma irradiation
- Bulk packed

PP = Polypropylene

Catalogue No	Alt. No	Material	Cap style	RCF, x g [max.]	Inner pack qty	Pack qty
CFT-900-011L	430897	PP	Plug seal	3,000	25	500
CFT-900-031F	430921	PP	Flat	3,000	25	500



## Tubes, centrifuge, 250mL/500mL

66 CORNING

### S

- Conical tubes for larger volume centrifugation
- Each case contains a recyclable foam rack to facilitate handling
- Max. RCF 6,000xg
- Sterile and certified non-pyrogenic
- Support cushion must be used unless centrifuge rotor has appropriate V-bottom

PP = Polypropylene

PEI = Polyethimide



Catalogue No	Alt. No	Description	Material	Cap style	RCF, x g [max.]	Inner pack qty	Pack qty
CFT-900-021Y	430776	250mL tube	PP	Plug	6,000	6	36
CFT-410-010A	431123	500mL tube	PP	Plug	6,000	6	36

### Accessory

Catalogue No	Alt. No	Description	Pack qty
CFT-896-502Q	430236	Support cushion for 250mL tube	6
CFT-410-030R	431124	Support cushion for 500mL tube	6

## Tubes, centrifuge, Thermo Scientific Nunc® EZ Flip™ cap

76

**Thermo**  
SCIENTIFIC

### S

- Integrated leakproof cap prevents cross contamination and lost closures
- Patented hinged-cap system
- Ergonomic design for one handed opening and closing
- Clarified polypropylene construction provides chemical resistance and high RCF tolerance
- White writing area and flat caps for sample identification
- Printed graduations
- Standard dimensions compatible with centrifuge rotors
- 15mL and 50mL versions are available

Thermo Scientific Nunc EZ Flip™ conical centrifuge tubes feature a patented, hinged flip-cap design in standard 15mL and 50mL sizes. These innovative tubes can be opened and closed with one hand, freeing the user's other hand for manipulations, such as pipetting.

Catalogue No	Alt. No	Description	RCF, x g [max.]	Inner pack qty	Pack qty
CFT-405-010A	362694	15mL bulk	8,500	50	500
CFT-405-030R	362695	15mL racked	8,500	50	500
CFT-405-050L	362696	50mL bulk	9,500	25	500
CFT-405-070F	362697	50mL racked	9,500	25	500

