



Tubes, microcentrifuge, polypropylene, conical

02

Kartell

Tubes for micro analysis. All 1.5mL tubes have attached hinged lids. Capping allows for safe transport, storage and centrifugation.

Catalogue No	Alt. No	Capacity, mL	Colour	Pack qty
TUL-645-030G	D0298-00	1.5	Neutral	1,000
TUL-645-050A	D0298-04	1.5	Blue	1,000
TUL-645-070R	D0298-06	1.5	Yellow	1,000
TUL-645-090L	D0298-12	1.5	Green	1,000
TUL-645-110Y	D0298-13	1.5	Pink	1,000
TUL-649-010L	D1299	0.5	Neutral	1,000



Tubes, microcentrifuge, polypropylene, Safe-Lock®

96

eppendorf



Designed for safe handling of valuable and toxic samples as well as radioactive substances or DNA.

- The small hook of the patented hinged lid clips around the rim of the tube giving trouble free heating
- Lid easily opened with one hand
- When open, tubes can be autoclaved at 121°C for 20 minutes
- Roughened writing surfaces (on side and lid)
- Biohazardous substances can be removed, free of aerosols, through a thin membrane in the centre of the lid with a cannula
- Autoclavable micropestle for resuspending pellets
- Available in Eppendorf® Biopur® guaranteed sterile, pyrogen-free, RNase-free, DNase-free, DNA-free, ATP-free
- Graduation scale
- Can be centrifuged at up to 23,500xg

These tubes for 0.5mL, 1.5mL and 2mL have a graduation scale and writing surfaces. The 0.5mL tube is optimised for use in PCR* experiments. Assorted colour packs are a mixture of yellow, red, blue, green and colourless.

Safe-Lock®

Catalogue No	Alt. No	Colour	Capacity, mL	Pack qty
TUL-150-010G	0030121023	Colourless	0.5	500
TUL-150-030A	0030121708	Assorted	0.5	500
TUL-150-050R	0030121112	Yellow	0.5	500
TUL-150-070L	0030121120	Red	0.5	500
TUL-150-090F	0030121139	Blue	0.5	500
TUL-150-110C	0030121147	Green	0.5	500
TUL-150-130T	0030121155	Amber	0.5	500
TUL-150-150N	0030120086	Colourless	1.5	1,000
TUL-150-170H	0030120694	Assorted	1.5	1,000
TUL-150-190B	0030120159	Yellow	1.5	1,000
TUL-150-210V	0030120167	Red	1.5	1,000
TUL-150-230P	0030120175	Blue	1.5	1,000
TUL-150-250J	0030120183	Green	1.5	1,000
TUL-150-270D	0030120191	Amber	1.5	1,000
TUL-150-290U	0030120094	Colourless	2.0	1,000
TUL-150-310R	0030120686	Assorted	2.0	1,000
TUL-150-330L	0030120205	Yellow	2.0	1,000
TUL-150-350F	0030120213	Red	2.0	1,000
TUL-150-370W	0030120221	Blue	2.0	1,000
TUL-150-390Q	0030120230	Green	2.0	1,000
TUL-150-410N	0030120248	Amber	2.0	1,000

PCR* clean Safe-Lock®

Catalogue No	Alt. No	Colour	Capacity, mL	Pack qty
TUL-155-030Q	0030123301	Natural	0.5	500
TUL-155-050K	0030123328	Natural	1.5	1,000
TUL-155-070E	0030123344	Natural	2.0	1,000

Biopur® Safe-Lock® individually sealed

Catalogue No	Alt. No	Colour	Capacity, mL	Pack qty
TUL-150-450B	0030121570	Natural	0.5	50
TUL-150-470S	0030121589	Natural	1.5	100
TUL-150-490M	0030121597	Natural	2.0	100

Accessory

Catalogue No	Alt. No	Description	Pack qty
TUL-150-430H	0030120973	Micropestle for 1.5mL and 2mL Safe-Lock® tubes and 1.7mL Safe-Twist screw cap tubes	10

*Polymerase Chain Reaction (PCR) is a process covered by patents owned by Hoffman-La Roche