

Tubes, microcentrifuge, polypropylene, snap cap, HydroLogix

03



**S**

- All HydroLogix tubes are certified RNase/DNase, RNA/DNA, ATP and pyrogen free
- Clear
- Graduation markings
- Frosted flat cap and side

Catalogue No	Alt. No	Capacity, mL	RCF, x g	Sterile, Yes/No	Pack qty
TUL-918-010X	3446	0.6	30,000	No	1,000
TUL-918-012K	3449	0.6	30,000	Yes	1,000
TUL-918-014G	3448	1.5	26,000	No	500
TUL-918-016C	3451	1.5	26,000	Yes	500
TUL-918-018V	3434	2.0	25,000	No	500
TUL-918-020L	3453	2.0	25,000	Yes	500

Tubes, microcentrifuge, polypropylene, locking lid, HydroLogix

03



**S**

- All HydroLogix tubes are certified RNase/DNase, RNA/DNA, ATP and pyrogen free
- Clear
- Graduation markings
- Frosted flat cap and side

These tubes feature an advanced locking system for safe and convenient centrifuging, boiling, freezing, incubation, or storage of samples for your most critical procedures.

The cap snaps closed, then the lock easily slips into place with one simple hand operation. The larger, rounder cap design provides more thumb contact while the smooth surface reduces the possibility of ripped or punctured gloves. TUL-918-034A, TUL-918-040F and TUL-918-046Q are individually wrapped.

Catalogue No	Alt. No	Capacity, mL	RCF, x g	Sterile, Yes/No	Pack qty
TUL-918-030Y	3454	0.6	30,000	No	1,000
TUL-918-032E	3455	0.6	30,000	Yes	1,000
TUL-918-034A	3455IW	0.6	30,000	Yes	500
TUL-918-036T	3456	1.5	24,000	No	500
TUL-918-038P	3457	1.5	24,000	Yes	500
TUL-918-040F	3457IW	1.5	24,000	Yes	500
TUL-918-042B	3458	2.0	25,000	No	500
TUL-918-044U	3459	2.0	25,000	Yes	500
TUL-918-046Q	3459IW	2.0	25,500	Yes	500

## Are you online?

- See the price you pay
- Live stock availability
- Order tracking
- Fully secure on line ordering

**Visit our website**

