



Tubes, microcentrifuge, polypropylene, screw-cap, conical, HydroLogix

03



Suitable for vapour phase liquid nitrogen storage, boiling and autoclaving.

- Made from a crystal clear resin, your samples are easily identifiable when using these tubes
- Available in amber for light sensitive applications
- Conical bottoms allow for more complete sample retrieval
- The superior O-ring design forms a positive seal against the rim of the screw cap to guard against leakage
- All HydroLogix tubes are certified RNase/DNase, RNA/DNA, ATP and pyrogen free
- Clear with frosted flat cap and side (TUL-918-058J, TUL-918-068G and TUL-918-078D are supplied without cap)
- Graduation markings

Catalogue No	Alt. No	Colour	Capacity, mL	RCF, x g	Sterile, Yes/No	Pack qty
TUL-918-050C	3460	Clear	0.5	25,000	No	500
TUL-918-052V	3460A	Amber	0.5	25,000	No	500
TUL-918-054R	3431	Clear	0.5	25,000	Yes	500
TUL-918-056N	3431A	Amber	0.5	25,000	Yes	500
TUL-918-058J	3462	Clear	0.5	25,000	No	500
TUL-918-060W	3464	Clear	1.5	20,000	No	500
TUL-918-062S	3464A	Amber	1.5	20,000	No	500
TUL-918-064X	3461	Clear	1.5	20,000	Yes	500
TUL-918-066K	3461A	Amber	1.5	20,000	Yes	500
TUL-918-068G	3466	Clear	1.5	25,000	No	500
TUL-918-070T	3468	Clear	2.0	20,000	No	500
TUL-918-072P	3468A	Amber	2.0	20,000	No	500
TUL-918-074L	3463	Clear	2.0	20,000	Yes	500
TUL-918-078D	3470	Clear	2.0	20,000	No	500



Tubes, microcentrifuge, polypropylene, screw-cap, free standing, HydroLogix

03



Suitable for vapour phase liquid nitrogen storage, boiling and autoclaving.

- Made from a crystal clear resin, your samples are easily identifiable when using these tubes
- Available in amber for light sensitive applications
- Conical bottoms allow for more complete sample retrieval
- The superior O-ring design forms a positive seal against the rim of the screw cap to guard against leakage
- All HydroLogix tubes are certified RNase/DNase, RNA/DNA, ATP and pyrogen free

With screw-cap and O-ring unattached

Catalogue No	Alt. No	Colour	Capacity, mL	RCF, x g	Sterile, Yes/No	Pack qty
TUL-918-080Q	3472	Clear	0.5	30,000	No	500
TUL-918-082M	3472A	Amber	0.5	30,000	No	500
TUL-918-090N	3474	Clear	1.5	18,000	No	500
TUL-918-092J	3474A	Amber	1.5	18,000	No	500
TUL-918-100N	3488	Clear	2.0	18,000	No	500
TUL-918-102J	3488A	Amber	2.0	18,000	No	500

With screw-cap and O-ring attached

Catalogue No	Alt. No	Colour	Capacity, mL	RCF, x g	Sterile, Yes/No	Pack qty
TUL-918-084Y	3465	Clear	0.5	30,000	Yes	500
TUL-918-086E	3465A	Amber	0.5	30,000	Yes	500
TUL-918-094F	3467	Clear	1.5	18,000	Yes	500
TUL-918-096B	3467A	Amber	1.5	18,000	Yes	500
TUL-918-104F	3469	Clear	2.0	18,000	Yes	500
TUL-918-106B	3469A	Amber	2.0	18,000	Yes	500

Without cap

Catalogue No	Alt. No	Capacity, mL	RCF, x g	Sterile	Pack qty
TUL-918-088A	3422	0.5	30,000	No	500
TUL-918-098U	3478	1.5	18,000	No	500



For Whatman VectaSpin™ filter centrifuge tubes please see Fisher Scientific Life Science Catalogue 2008-2009, page 676