

Tube racks, microcentrifuge, polypropylene, 4-way

Fisherbrand®

41

a



Four way heavy duty racks feature a unique system of tabs and slots which when linked together create a large rack.

Each rack can hold 4 x 50mL conical tubes, 12 x 15mL conical tubes, 32 x 1.5mL microtubes or 32 x 0.5mL microtubes. Rack measures 174mm x 95mm x 52mm. Autoclavable.


Catalogue No	Description	Pack qty
FB70960	Microtube rack 4-way, natural	5
FB70961	Microtube rack 4-way, blue	5
FB70962	Microtube rack 4-way, green	5
FB70963	Microtube rack 4-way, pink	5
FB70964	Microtube rack 4-way, yellow	5
FB70965	Microtube rack 4-way, orange	5
FB70966	Microtube rack 4-way, assorted colours (blue, green, pink, yellow and orange)	5

Tube racks, microcentrifuge, polypropylene, 20 well

Fisherbrand®

41

a



Polypropylene racks for storing 20 tubes. Racks lock together on each side.

Tubes are easily removed using 'hold-it' finger slot on the front of the rack. Racks measure 213mm x 90mm x 50mm. Autoclavable.

Catalogue No	Description	Pack qty
FB70971	Microtube rack 20 well, blue	6
FB70973	Microtube rack 20 well, white	6

Tube rack, microcentrifuge, polypropylene, 3-way

Fisherbrand®

41

a



Three way polypropylene racks connect to one another on either end. Convenient finger slots are moulded onto the rack and each of the three tiers accommodates a different size tube. Lowest level has 24 wells to hold individual 0.2mL tubes or strips of 8 or 12 tubes. The middle tier has 14 wells to hold 0.5mL tubes. The top tier has 12 wells to hold 1.5mL tubes. Autoclavable. Rack measures 164mm x 123mm x 62mm and is supplied in assorted colours (blue, green, purple, yellow, and orange).

Catalogue No	Description	Pack qty
FB70969	Microtube rack 3-way, assorted colours	6

Tube racks, microcentrifuge, polypropylene, cube

Fisherbrand®

41

a



User friendly cube racks.

Slide the polypropylene racks together in multiple orientations - work with different sized tubes at the same time. Each rack can hold 4 x 50mL tubes, 10 x 15mL conical tubes, 12 x 75mL or 100mL tubes or 16 x 1.5mL to 2.0mL tubes. Autoclavable. Racks measure 108mm x 113mm x 113mm.

Catalogue No	Description	Pack qty
FB70975	Cube rack, assorted colours (one each of blue, green, yellow, orange, and pink)	5
FB70976	Cube rack, natural	5