


Centrifuge tubes racks

Tube racks, microcentrifuge, polypropylene, 80 well

41

Fisherbrand®

a



Polypropylene 80 well microtube racks ideal for freezer storage.

Racks hold 1.5mL to 2mL microtubes in a 5 x 16 array. Tubes sit in large wells with conical bottoms. Racks measure 225mm x 67mm x 28mm. Autoclavable.

Catalogue No	Description	Pack qty
FB71000	Natural	5
FB71002	Blue	5
FB71004	Green	5
FB71006	Pink	5
FB71008	Yellow	5
FB71010	Orange	5
FB71012	Assorted colours (one of each colour)	5

Tube racks, microcentrifuge, polypropylene, 100 well

41

Fisherbrand®

a



Durable 100 well storage racks feature three point hinges and clasp for secure lid closure.

Hold 1.5mL to 2mL microtubes and cryovials. User friendly racks feature moulded grid lines on the lid, imprinted coordinates on the bottom of the base and on the box. Racks measure 141mm x 151mm x 57mm. Autoclavable.

Catalogue No	Description	Pack qty
FB71020	100 well rack, natural	5
FB71021	100 well rack, assorted colours (blue, green, purple, yellow and orange)	5

Tube racks, centrifuge, Rota-Rack™

41

Fisherbrand®

a



Each module can be independently rotated enabling the use of several sizes of test tubes, culture tubes, centrifuge tubes or microcentrifuge tubes at the same time. Click each module in to place, load with tubes and carry to a new location. The racks remain fully orientated until they are twisted again.

Each module of the large Rota-Rack™ holds 2 x 50mL tubes, 8 x 15mL tubes, 6 x 16mm to 20mm tubes, or 10 x 10mm to 13mm tubes.

Each module of the small Rota-Rack™ holds 6 x 15mL tubes, 9 x 1.5/2mL tubes, 12 x 0.5/0.6mL tubes or 32 independent 0.2mL PCR* tubes or 4 x 8 tube strips.

Each rack has modules in green, pink, blue and yellow, and is fully autoclavable.

Catalogue No	Description	Dimensions (w x d x h), mm
FB51101	Rota-Rack™, large	224 x 112 x 112
FB51102	Rota-Rack™, small	208 x 100 x 100

Tube racks, microcentrifuge, polypropylene, 50 well

41

Fisherbrand®

a



Rugged polypropylene racks have moulded grid lines on the hinged lid, imprinted coordinates at the base of each well, allowing easy identification.

Stackable 50 well racks hold 1.5mL to 2mL microtubes and cryovials. Racks measure 141mm x 92mm x 56mm. Autoclavable.

Catalogue No	Description	Pack qty
FB71023	50 well rack, natural	5
FB71024	50 well rack, assorted colours (blue, green, purple, yellow and orange)	5