

Tube racks, microcentrifuge, polypropylene, PCR\* tube, reversible

**Fisherbrand®**

41

**a**



Reversible polypropylene rack with removable hinged lid that snaps in place on either side of the rack.

PCR\* side of the rack has 168 wells that hold 0.2mL tubes. The wells accommodate 8 or 12 tube strips or individual tubes. The opposite side of the rack has 40 wells that hold 0.5mL tubes. There are 12 wells for 1.5mL tubes on both sides of the rack. Autoclavable. Rack measures 206mm x 131mm x 54mm. Supplied in a pack of assorted colours.

Catalogue No	Description	Pack qty
FB71030	Reversible PCR* rack, assorted colours (blue, green, purple, yellow and orange)	5

Tube racks, microcentrifuge, polypropylene, PCR\* tube

**Fisherbrand®**

41

**a**



Sturdy polypropylene racks ideal for pre- or post-PCR\* sampling applications.

With removable hinged lid accommodates individual 0.2mL tubes or strips of 8 or 12 tubes.

Tube wells are easily identified with imprinted numbers and letters. Racks measure 130mm x 98mm x 33mm. Autoclavable.

Catalogue No	Description	Pack qty
FB71035	96 well PCR* rack, assorted colours (blue, green, pink, yellow and orange)	5
FB71037	96 well PCR* rack, natural	5
FB71039	96 well PCR* rack, blue	5
FB71041	96 well PCR* rack, green	5
FB71043	96 well PCR* rack, yellow	5

Tube rack, microcentrifuge, polypropylene, PCR\* transfer tray

**Fisherbrand®**

41

**a**



Save time when working with individual 0.2mL tubes or strips of 8 or 12 tubes.

Transfer tray is a thin strip of plastic designed to rest perfectly over a 96 well tube rack. Just use handles to lift and carry, ideal when using with thermal cyclers. Autoclavable. Tray measures 114mm x 74mm x 11mm. Supplied in a pack of assorted colours.

Catalogue No	Description	Pack qty
FB71050	PCR* transfer tray, assorted colours (blue, green, purple, yellow and orange)	5

Tube racks, microcentrifuge, polypropylene, 96 well, reversible

**Fisherbrand®**

41

**a**



Polypropylene racks have 96 wells to accommodate 0.5mL microtubes on one side, or 96 wells for 1.5mL to 2mL microtubes on the other side.

Arranged in standard 8 x 12 array. Handy lid snaps into place in the tabs on each end of rack. Grid lines are moulded on the side of the rack making it easier to see the wells. Raised rings are moulded on the opposite side of the rack making it easier to replace the 1.5mL to 2mL tubes from these wells. Autoclavable. Racks measure 246mm x 121mm x 50mm. Available in natural or assorted colours.

Catalogue No	Description	Pack qty
FB71055	96 well reversible rack, assorted colours (blue, green, pink, yellow and orange)	5
FB71057	96 well reversible rack, natural	5
FB71059	96 well reversible rack, blue	5
FB71061	96 well reversible rack, green	5
FB71063	96 well reversible rack, yellow	5