

Vials, Shell (neckless vials)

17



For use in all non-robotic autosamplers where the vial sits in the carousel.

Suitable for Waters 96, 8mm x 36mm

Catalogue No	Alt. No	Description	Capacity, mL	Pack qty
VGA-350-010A	1-NWV-C	Clear with polyethylene cap	1	200
VGA-350-020U	1-NWV(A)-C	Amber with polyethylene cap	1	200
VGA-355-010Q	1-NWV	Clear without cap	1	500

Plugs, polyethylene

Catalogue No	Alt. No	Description	Pack qty
VGA-358-020S	8-NPWP	8mm, for Waters 96 position shell vials	1,000

Suitable for Waters 48, 15mm x 46mm

Catalogue No	Alt. No	Description	Capacity, mL	Pack qty
VGA-350-030R	4-NWV-C	Clear with polyethylene cap	4	100



Inserts and insert supports

17



Use 02-MTV with 2-SV, 2-SV(A) and 2SVG, 2-SVV, 2SVW(A) vials.

Use 02-MTVGW with 2-CV, 2-CV(A), 2-CVG, 2-DV, 2-RV vials.

Catalogue No	Alt. No	Description	Capacity, µL	Pack qty
⑤ VGA-100-500L	02-MTV	Clear, standard	200	1,000
VGA-100-504D	02-MTVGW	Clear, wide neck	200	1,000
⑥ VGA-100-502L	MTS-1	Polyethylene, support legs, self-centering for inserts	-	500



Vials and caps, headspace, crimp top, clear, 20mm

17



- Vials and caps sold separately

Technical Specification - Specific

	VGA-100-050Q	VGA-100-070K	VGA-100-080H	VGA-100-085U	VGA-100-090E
Capacity, mL	6	10	20	22	27
Instrument	Perkin Elmer	Varian, Tekmar	Agilent, Dani, Dani Purge & Trap, Fisons HS 850, Perkin Elmer	Variant, Tekmar, Al Cambridge	Shimadzu HSS-2B

Clear glass vials

Catalogue No	Alt. No	Dimensions, mm	Capacity, mL	Cap size, mm	Pack qty
VGA-100-050Q	6-CV	22 x 38	6	20	125
VGA-100-070K	10-CV	22 x 45	10	20	125
VGA-100-080H	20-CV	22 x 75	20	20	125
VGA-100-085U	22-CV	22 x 75	22	20	125
⑦ VGA-100-090E	27-CV	30 x 56	27	20	100

Blank caps, without seal

Catalogue No	Alt. No	Description	Dimensions, mm	Pack qty
VGA-104-900R	20-ACB	Aluminium	20 x 7	500
VGA-104-902N	20-MCB	Magnetic, tin plate	20 x 7	500
VGA-104-906F	20-MCBC	Composite blue, magnetic	20 x 7	500

Seals

Catalogue No	Alt. No	Description	Dimensions, mm	Pack qty
VGA-104-910X	20-CB3	Chlorobutyl	20 x 3	1,000
VGA-104-924D	20-ST3	Silicone/PTFE	20 x 3	500

Aluminium caps with seals, pre-assembled

Catalogue No	Alt. No	Description	Dimensions, mm	Pack qty
VGA-105-010V	20-AC-CBT3	Chlorobutyl/PTFE seal	20 x 7	500
VGA-105-020S	20-AC-ST3	Silicone/PTFE seal	20 x 7	500

Plugs

Catalogue No	Alt. No	Description	Dimensions, mm	Pack qty
VGA-104-921J	20-B3P	Butyl rubber plug	20 x 3	500

