



Vials, shell, 1mL, 2mL and 4mL with polyethylene plug

17

Fisherbrand®

NEW

Clear and amber Fisherbrand shell vials are suitable for use with Alcott, Gilson, Shimadzu and Waters® (Wisp 96 and 48 position carousel) autosamplers.

Clear glass with transparent polyethylene plug

Catalogue No	Capacity, mL	Dimensions, mm	Pack qty
FB67171	1.0	7.8 x 35	100
FB67173	1.0	8 x 40	100
FB67177	2.0	12 x 32	100
FB67181	4.0	15 x 45	100

Amber glass with transparent polyethylene plug

Catalogue No	Capacity, mL	Dimensions, mm	Pack qty
FB67175	1.0	8 x 40	100
FB67179	2.0	12 x 32	100
FB67183	4.0	15 x 45	100



Vials, crimp top and micro, 8mm

17

Fisherbrand®

NEW


Fisherbrand clear and amber glass crimp top and microvials are suitable for use with the following autosamplers: Agilent; Beckman; Carlo Erba; CTC; Fisons; Perkin Elmer; Shimadzu; Thermo and VWR (Merck®)/Hitachi, etc.

Clear glass

Catalogue No	Description	Capacity, mL	Dimensions, mm	Pack qty
FB67193	Round bottom	0.3	5.5 x 31.5	100
FB67197	Conical bottom	0.3	5.5 x 31.5	100
FB67207	Conical bottom	0.7	7 x 40	100
FB67199	Flat bottom	0.7	7 x 40	100
FB67191	Flat bottom	0.8	8.2 x 30	100
FB67203	Flat bottom	1.2	8.2 x 30	100

Amber glass

Catalogue No	Description	Capacity, mL	Dimensions, mm	Pack qty
FB67195	Conical bottom	0.7	7 x 30	100
FB67209	Conical bottom	0.7	7 x 40	100
FB67201	Flat bottom	0.7	7 x 40	100
FB67205	Flat bottom	1.2	8.2 x 40	100



Crimp seals, aluminium, 8mm

17

Fisherbrand®

NEW

Catalogue No	Description	Thickness, mm	Pack qty
FB67211	With virgin PTFE seal	0.25	100
FB67213	With red-orange natural rubber/transparent TEF seal	1.0	100
FB67215	With red-orange natural rubber/transparent TEF seal (Agilent quality)	1.0	100
FB67217	With red rubber, beige seal (Agilent quality)	1.0	100
FB67221	With red PTFE/white silicone/red PTFE seal	1.0	100
FB67219	With white silicone/red PTFE seal	1.3	100
FB67223	With white silicone/red PTFE seal with slit	1.3	100



Inserts, microvial, for crimp vials, 11mm (wide opening)

17

Fisherbrand®

NEW

Catalogue No	Description	Capacity, mL	Dimensions, mm	Pack qty
FB67231	Clear glass, 15mm top	0.1	6 x 31	1,000
FB67233	Clear glass, 12mm top	0.1	6 x 31	1,000
FB67235	Clear glass, with assembled plastic spring	0.1	5.7 x 29	1,000
FB67237	Clear glass, flat bottom	0.2	6 x 31	1,000