

Vials, crimp top and micro, 2mL, 11mm (wide opening)

Fisherbrand

NEW

17

Fisherbrand amber and clear glass crimp top and micro vials, suitable for use with the following autosamplers: Agilent, Carlo Erba, CTC, Dani, Fisons, Gerstel, Jasco, Perkin Elmer, Shimadzu, Spark, Thermo, Varian, etc.

Clear glass

Catalogue No	Description	Capacity, mL	Dimensions, mm	Pack qty
FB67249	Crimp vial with integrated 0.2mL micro insert, 1st hydrolytic class, graduated with marking spot	0.2	12 x 32	100
FB67253	TopSert™TPX snap ring vial with integrated 0.2mL micro insert	0.2	12 x 32	100
FB67247	Microvial, 15µL, conical	1.1	12 x 32	100
FB67251	Microvial, conical bottom	1.1	12 x 32	100
FB67241	Wide opening	2	12 x 32	100
FB67243	Wide opening, graduated with marking spot	2	12 x 32	100

Amber glass

Catalogue No	Description	Capacity, mL	Dimensions, mm	Pack qty
FB67245	Amber glass, wide opening, graduated with marking spot	2	12 x 32	100



Caps, crimp, aluminium, 11mm

Fisherbrand

NEW

17

Catalogue No	Description	Thickness, mm	Pack qty
FB67261	With red-orange natural rubber/transparent TEF seal (Agilent quality)	1.0	100
FB67265	With red rubber/PTFE seal (Agilent quality)	1.0	100
FB67269	With red PTFE/white silicone/red PTFE seal	1.0	100
FB67263	With red-orange natural rubber/transparent TEF seal	1.3	100
FB67267	With white silicone/red PTFE seal	1.3	100
FB67271	With white silicone/blue PTFE seal, with cross-slit	1.5	100



Caps, crimp, magnetic, 11mm


Fisherbrand

NEW

17

- For CTC Pal+TriPlus autosampler

Catalogue No	Description	Thickness, mm	Pack qty
FB67277	With red PTFE/white silicone/red PTFE seal	1.0	100
FB67275	With white silicone/red PTFE seal	1.3	100



Vials, snap ring, 2mL and microvials, 11mm (wide opening)

Fisherbrand

NEW

17

Clear and amber glass Fisherbrand snap ring and microvials are suitable for use with the following autosamplers: Agilent, CTC, Dani, Dionex, Jasco, Shimadzu, Spark, Thermo, Varian, VWR (Merck®)/Hitachi, Waters®, etc.

Clear glass

Catalogue No	Description	Capacity, mL	Dimensions, mm	Pack qty
FB67327	TopSert™PX snap ring vial with integrated 0.2mL glass micro insert	0.2	12 x 32	100
FB67321	Clear glass, wide opening	2	12 x 32	100
FB67323	Clear glass, wide opening, graduated with marking spot	2	12 x 32	100

Amber glass

Catalogue No	Description	Capacity, mL	Dimensions, mm	Pack qty
FB67325	Amber glass, wide opening, graduated with marking spot	2	12 x 32	100

