



Micro inserts, for 2mL snap ring vials (wide opening)

Fisherbrand

17

NEW

Catalogue No	Description	Capacity, mL	Dimensions, mm	Pack qty
FB67311	Clear glass, 15mm top	0.1	6 x 31	1,000
FB67315	Clear glass, with assembled plastic spring	0.1	5.7 x 29	1,000
FB67313	Clear glass, 12mm top	0.2	7 x 31	1,000
FB67317	Clear glass, flat bottom	0.2	6 x 31	1,000



Seals, snap ring, polyethylene, 11mm with 6mm centre hole, ND11

Fisherbrand

17

NEW

Further ND11 snap ring seals, hard cap version, with transparent, blue, red, green and yellow caps are available on request.

Transparent cap, soft cap version

Catalogue No	Description	Thickness, mm	Pack qty
FB67281	With red-orange natural rubber/transparent TEF seal	1.0	100
FB67283	With red rubber/PTFE seal (Agilent quality)	1.0	100
FB67287	With red PTFE/white silicone/red PTFE seal	1.0	100
FB67289	With white silicone/blue PTFE seal, with cross-slit	1.0	100
FB67291	With white silicone/red PTFE seal, pre-cut (Y)	1.0	100
FB67285	With white silicone/red PTFE seal	1.3	100

Blue cap, soft cap version

Catalogue No	Description	Thickness, mm	Pack qty
FB67293	With red-orange natural rubber/transparent TEF seal	1.0	100
FB67295	With red rubber/PTFE seal (Agilent quality)	1.0	100
FB67299	With red PTFE/white silicone/red PTFE seal	1.0	100
FB67301	With white silicone/blue PTFE seal, with cross-slit	1.0	100
FB67303	With white silicone/red PTFE seal, pre-cut (Y)	1.0	100
FB67297	With white silicone/red PTFE seal	1.3	100



Vials, screw and micro, 8mm (small opening)

17

Fisherbrand

NEW

Fisherbrand screw top vials and microvials are for use with the following autosamplers: Beckman, CTC, Gilson, Knauer, Shimadzu, Spark, Varian, VWR (Merck®)/Hitachi, etc.

Clear glass

Catalogue No	Description	Capacity, mL	Dimensions, mm	Pack qty
FB67351	Vial, small opening	2	12 x 32	100
FB67353	Vial, small opening, graduated with marking spot	2	12 x 32	100
FB67357	Microvial, small opening	2	12 x 32	100

Amber glass

Catalogue No	Description	Capacity, mL	Dimensions, mm	Pack qty
FB67355	Small opening, graduated with marking spot	2	12 x 32	100



Micro inserts, for 2mL screw vials (small opening)

17

Fisherbrand

NEW

Catalogue No	Description	Capacity, mL	Dimensions, mm	Pack qty
FB67361	Clear glass, 15mm top	0.1	5 x 31	1,000
FB67363	Clear glass, 9mm top	0.1	5 x 31	1,000
FB67365	Clear glass, with assembled plastic spring	0.1	5 x 29	1,000
FB67367	Clear glass, flat bottom	0.2	5 x 31	1,000