



Vials, EPA, screw, 24mm opening

17

Fisherbrand

NEW

EPA vials are for use with Agilent, Dionex, Shimadzu, Tekmar, Thermo, Varian, etc. autosamplers.

Clear glass

Catalogue No	Capacity, mL	Dimensions, mm	Pack qty
FB67491	20	27.5 x 57	100
FB67495	30	27.5 x 72.5	100
FB67499	40	27.5 x 95	100
FB67503	60	27.5 x 140	100

Amber glass

Catalogue No	Capacity, mL	Dimensions, mm	Pack qty
FB67493	20	27.5 x 57	100
FB67497	30	27.5 x 72.5	100
FB67501	40	27.5 x 95	100
FB67505	60	27.5 x 140	100



Caps, screw, polypropylene, with septa, 24mm, centre hole and closed top, UltraBond™

17

Fisherbrand

NEW

UltraBond™ caps and septa form an inseparable unit so that the septum liner cannot fall out or be pushed into the vial when being penetrated with a blunt needle. The 100% firm fit is achieved by a patented process without the use of any glue or adhesive, but by changing the molecular structure of the surface of the cap and septa.

Catalogue No	Description	Septa material	Thickness, mm	Pack qty
FB67507	White screw cap, 15mm centre hole	Natural silicone/beige PTFE (EPA quality)	3.2	100
FB67509	White screw cap, closed top	Natural silicone/beige PTFE (EPA quality)	3.2	100



Vials, headspace, precision thread, 18mm opening

17

Fisherbrand

NEW

Precision thread headspace vials are compatible with CTC PAL, Varian, Gerstel, Atas, Shimadzu and Agilent autosamplers.

Clear glass

Catalogue No	Description	Capacity, mL	Dimensions, mm	Pack qty
FB67511	Precision thread headspace vial, clear glass, round bottom	10	22.5 x 46	100
FB67515	Precision thread headspace vial, clear glass, round bottom	20	22.5 x 75.5	100

Amber glass

Catalogue No	Description	Capacity, mL	Dimensions, mm	Pack qty
FB67513	Precision thread headspace vial, amber glass, round bottom	10	22.5 x 46	100
FB67517	Precision thread headspace vial, amber glass, round bottom	20	22.5 x 75.5	100