

HPLC columns, Thermo Scientific Hypersil™ BDS

16



- Base deactivated and end capped for reduced peak tailing and excellent peak symmetry
- Highly reproducible
- Robust columns with long lifetimes
- Excellent performance with basic, neutral and acidic compounds

C₁₈

Catalogue No	Alt. No	Column I.D., mm	Particle size, µm	Length, mm
HSA-117-505L	28103-054630	4.6	3	50
HSA-217-510K	28105-104630	4.6	5	100
HSA-217-515A	28105-154630	4.6	5	150
HSA-217-525U	28105-254630	4.6	5	250

C₁₈ guard columns

Catalogue No	Alt. No	Column I.D., mm	Particle size, µm	Length, mm	Pack qty
HSA-117-600R	28103-014001	4.0 and 4.6	3	10	4
HSA-217-600J	28105-014001	4.0 and 4.6	5	10	4

C₈

Catalogue No	Alt. No	Column I.D., mm	Particle size, µm	Length, mm
HSA-123-505F	28203-053030	3.0	3	50
HSA-223-510E	28205-104630	4.6	5	100
HSA-223-515R	28205-154630	4.6	5	150
HSA-223-525X	28205-254630	4.6	5	250

C₈ guard columns

Catalogue No	Alt. No	Column I.D., mm	Particle size, µm	Length, mm	Pack qty
HSA-123-600L	28203-013001	3.0	3	10	4
HSA-223-600D	28205-014001	4.0 and 4.6	5	10	4

Cyano CN

Catalogue No	Alt. No	Column I.D., mm	Particle size, µm	Length, mm
HSA-248-515T	28805-154630	4.6	5	150
HSA-248-525Q	28805-254630	4.6	5	250

Cyano CN guard column

Catalogue No	Alt. No	Column I.D., mm	Particle size, µm	Length, mm	Pack qty
HSA-248-600F	28805-014001	4.0 and 4.6	5	10	4

Phenyl

Catalogue No	Alt. No	Column I.D., mm	Particle size, µm	Length, mm
HSA-243-515D	28905-154630	4.6	5	150
HSA-243-525A	28905-254630	4.6	5	250

Phenyl guard column

Catalogue No	Alt. No	Column I.D., mm	Particle size, µm	Length, mm	Pack qty
HSA-243-600M	28905-014001	4.0 and 4.6	5	10	4

Silicas

A wide range of silicas is available to suit differing requirements from the bench through to product purification. The range of particle sizes includes 35-70µm and 40-63µm, 60Å and 100Å pore sizes.

See the Fisher Scientific Chemicals catalogue for the full range of products.

DAVISIL