



## Flash systems, accessories and spares

- Highly resistant polypropylene pre-packed cartridges
- Leak-free connection system, even at pressures up to 10bar
- Supplied with Lot Number and Quality Certificate
- All cartridges are individually sealed
- Supplied in robust shock absorbent case

The following accessories are the most commonly requested Sepacore® items. Details of self-pack cartridges available, please contact Customer Service. Please note that cartridges must not be used with pressures higher than 10bar.

### 1 Pre-packed cartridges

Catalogue No	Alt. No	Diameter, mm	Silica	Pack qty
CKH-675-050L	044884	12	4g	20
CKH-675-060Y	044885	12	8g	20
CKH-675-070F	044886	40	40g	10
CKH-675-080C	044887	40	90g	10
CKH-675-200S	054860	12	4g	10
CKH-675-210P	054861	12	8g	10
CKH-675-220M	054862	40	40g	5
CKH-675-230J	054863	40	90g	5



### Sample loading accessories

Catalogue No	Alt. No	Description	Pack qty
CKH-601-500K	044850	Injection valve 1/4 inch UNF-28 for injection of sample volumes up to 5mL	1
CKH-601-510H	044851	6 way valve. Enables range of other sample loading accessories to be used	1
CKH-601-520E	044852	Sample loop for sample volume 1mL to 20mL, requires 6 way injection unit	1
CKH-690-030K	44030	PrepElut glass column for loading dry or pre-adsorbed samples, 18mL capacity, requires 6 way injection unit	1
CKH-690-075L	044877	PrepElut adapter for loading dry/pre-adsorbed samples up to 8mL or 17mL onto 12mm cartridges, requires 6 way valve	1



### Collection accessories

Catalogue No	Alt. No	Description	Pack qty
CKH-660-512N	54293	Column module	1
2 CKH-660-500U	44957	Rack FC12 - for 12 x 40mm dia. tubes (tubes not included)	4
CKH-660-530L	44960	Glass tubes, 250mL for FC12	10
CKH-660-510R	44958	Rack FC30 - for 30 x 26mm dia. tubes (tubes not included)	4
CKH-660-540Y	44961	Glass tubes, 50mL for FC30	20
CKH-660-520X	44959	Collection rack FC60 - for 60 x 17mm diameter tubes (tubes not included)	4
CKH-660-550F	44962	Glass tubes, 16mm diameter, 160mm height, 20mL capacity for FC60 collection rack	100
CKH-601-564H	45288	Waste level sensor (float type) for CKH-615-010T. GL45 thread. Detection range 0mL to 10mL	1
CKH-601-565F	45285	Solvent level sensor, float type. Detects how full/empty solvent/waste reservoir bottle is. Not for use with corrosive solvents such as strong acids, strong bases or concentrated salt solutions	1
3 CKH-601-620A	44021	Tubing fitting kit	1
CKH-601-585W	54400	Tool for easy removal of grey or green fitting	1
CKH-601-690C	44340	Solvent filter. Prevents dust and particulates from reaching the pump	1
CKH-999-010H	16422	SepacoreRecord standard licence - includes SepacoreRecord data acquisition and SepacoreRecord fraction collector licence	1