

Disinfectant wipes, Azowipes

50



Bactericidal hard surface disinfectant wipes.

200 wipes per tub.
Use biocides safely. Always read the label and product information before use.



Catalogue No	Dimensions, mm
HYG-231-001T	200 x 200

Wipes, cleanroom, Class 10, ISO Class 4 Super Polx® 1200

04




- Continuous filament, non-run knitted polyester
- Laundered and packaged in a Class 10 cleanroom
- Extremely low in particle and fibre generation
- Low extractables
- Manufactured and tested with strict adherence to quality assurance guidelines
- Vacuum sealed

Catalogue No	Dimensions, mm	Pack qty
CMH-600-090K	230 x 230	150 wipes
CMH-600-100K	300 x 300	75 wipes

Cloths

03



	Catalogue No	Description	Pack qty
1	CMC-482-M	Wiping cloth, resin bonded cellulose fibre, soft in texture and absorbent, may be used dry or moistened for wiping and cleaning purposes, will withstand careful washing, (356mm x 343mm)	25
2	CMC-450-W	Cleaning pads, very fine grade for cleaning, scouring and sundry uses, (240mm x 150mm)	10
3	CMC-300-T	Drying linen/union, long lasting, absorbent and strong, (760mm x 480mm)	3
4	CMC-330-V	General purpose, for general cleaning use, 100% cotton, assorted colours, (450mm x 700mm)	3
	CMC-360-A	Polishing cloth, high quality, all cotton, soft, close texture, ideal for use with commercial polishes (500mm x 500mm)	5
	CMC-510-T	Swabbing cloth, for wiping down laboratory surfaces and general spillages, open mesh, (500mm x 450mm)	10

Contamination control mat, multi-layer, disposable, white

04



Low profile. Pack of 8 tack mats x 30 layers.

Catalogue No	Dimensions, mm	Pack qty
CMK-940-020S	457 x 1,168 (approx.)	8

Face masks and dispenser, cleanroom, Class 100 ISO Class 5

04



Disposable masks specifically designed for Class 100 cleanrooms.

Very low particulate generation; filtration efficiency at 1.0mm is 96%, consisting of three layers of lint-free material. The two outer layers are made from lint-free, wet laid, non-woven material and the centre filtration medium is made of fine, melt blown polypropylene. Each mask has a metal nose piece which is easy to bend providing a contoured comfortable fit. Folded portion of the mask extends to cover chin area. Individually packed.

The dispenser can be usefully stationed in the changing area of any cleanroom.

Catalogue No	Size	Pack qty
CMK-530-010X	Standard	50
CMK-530-020L	Wide	50

Accessory

Catalogue No	Description
CMK-530-500U	Dispenser for CMK-530-010X, holds 50 masks