

Containers, polystyrene, bijou, with coverslip, 7mL 56



Ideal for work on the culture of intracellular organisms (e.g. *Chlamydia trachomatis*, *Rickettsia* etc.).

- Coverslip enables the growth of a monolayer of cells into which the intracellular organism has been infected
- Plastic cap
- Coverslip can be removed for examination/study
- Aseptically manufactured

Catalogue No	Alt. No	Description	Pack qty
DIS-080-170S	129AX/1	No label	700
DIS-085-170L	129BX/1	Plain label	700



Containers, polystyrene, 7 to 250mL 56



A range of high quality polystyrene sample containers, manufactured under strict aseptic conditions. All Sterilin polystyrene containers are CE marked to comply with the In Vitro Diagnostic Directive. Aseptically manufactured from cytotoxicity tested, FDA approved virgin polystyrene. Containers are leak tested to EN 14254 Annex D. Most sizes have a choice of metal or plastic caps.

7mL bijou

Catalogue No	Alt. No	Label type	Cap material	Pack qty
DIS-080-010R	129A	No	PE	700
DIS-085-010K	129B	Plain	PE	700

30mL universal

Catalogue No	Alt. No	Description	Label type	Cap material	Pack qty
DIS-080-100Q	128A	FS cap	None	PP	400
DIS-080-110N	128SA	Spoon	None	PE	400
DIS-080-120K	128B-BAC	Boric Acid FS cap	Yes	PP	400
DIS-085-022D	128SC	Spoon	Plain	PE	400
DIS-085-100J	128B	FS cap	Printed	PP	400
DIS-085-105W	128C	FS cap	Plain	PP	400
DIS-085-110G	128SB	Spoon	Printed	PE	400

60mL containers

Catalogue No	Alt. No	Description	Label type	Cap material	Pack qty
DIS-080-024G	125AM	FS cap	None	M	300
DIS-080-025E	125AP	Plastic cap	None	PE	300
DIS-085-024W	125BM	FS cap	Printed	M	300
DIS-085-025U	125BP	Plastic cap	Printed	PE	300
DIS-085-028X	125CP	Plastic cap	Plain	PE	300
DIS-085-120D	125CM	FS cap	Plain	M	300

100mL containers

Catalogue No	Alt. No	Description	Label type	Cap material	Pack qty
DIS-080-030L	185AM	FS cap	None	M	200
DIS-080-031J	185AP	Plastic cap	None	PE	200
DIS-085-030E	185BM	FS cap	Printed	M	200
DIS-085-031C	185BP	Plastic cap	Printed	PE	200
DIS-085-035R	185CM	FS cap	Plain	M	200
DIS-085-036P	185CP	Plastic cap	Plain	PE	200

150mL containers

Catalogue No	Alt. No	Description	Label type	Cap material	Pack qty
DIS-080-040Y	165A	FS cap	None	M	120
DIS-085-040B	165B	FS cap	Printed	M	120
DIS-085-050V	165C	FS cap	Plain	M	120

