

Vials, cryogenic, internal and external thread, CryoTube™, Thermo Scientific Nunc®



76

Thermo
SCIENTIFIC


Designed for cryogenic storage of cell cultures, microbiological cultures, sperm, blood and serum specimens.

- Perform well in mechanical freezers, and in liquid nitrogen (gas and liquid¹ phase)
- Internal thread with a silicone gasket provides the best possible seal
- For extreme storage conditions like liquid nitrogen CryoFlex™ must be used¹

¹WARNING: Do not use CryoTubes for freezing in liquid nitrogen (liquid phase storage)



▲ Stoppers with silicon gasket for tubes with internal thread

● Screw cap for tubes with external thread

*Closure design does not permit the use of CryoColourcode



Internal thread, free standing, 12.5mm (max.) diameter, polypropylene tubes and screw caps

Catalogue No	Alt. No	Bottom	Volume, working, mL	Starfoot [Y/ N]	Length, mm	Diameter, mm	Inner pack qty	Pack qty
⑨ CRY-960-010T	363401	U	1.8	No	48	12.5 (max.)	50	500
① CRY-960-030N	366524	U	3.6	No	70	12.5 (max.)	100	400
② CRY-960-050H	363452	U	4.5	No	92	12.5 (max.)	50	300
⑤ CRY-960-070B	366656	V	1.0	No	42	12.5 (max.)	50	500
④ CRY-960-090S	368632	U	1.8	No	49	12.5 (max.)	50	450

Internal thread, starfoot design, 12.5 (max.) diamter, polypropylene tubes and screw caps

Catalogue No	Alt. No	Bottom	Volume, working, mL	Starfoot [Y/ N]	Length, mm	Inner pack qty	Pack qty
⑥ CRY-960-110P	377224	V	1.0	Yes	42	50	500
⑧ CRY-960-130J	377267	U	1.8	Yes	49	50	450
③ CRY-960-150D	379189	U	3.6	Yes	72	100	400
⑦ CRY-960-170U	379146	U	4.5	Yes	92	300	300

External thread, starfoot design, 12.5mm diameter, polypropylene tubes, polyethylene screw caps

Catalogue No	Alt. No	Bottom	Volume, working, mL	Starfoot [Y/ N]	Inner pack qty	Pack qty
⑩ CRY-960-190X	375353	V	1.0	Yes	41	500
⑪ CRY-960-210L	375418	U	1.8	Yes	48	450
⑫ CRY-960-230F	337516	U	4.5	Yes	91	300
⑬ CRY-960-250W	347597	V	1.0	No	41	500
⑭ CRY-960-251U	347627	U	1.8	No	48	450
⑮ CRY-960-253Q	347643	U	4.5	No	91	300

External thread 12.4mm diameter, polypropylene tubes, polyethylene screw caps* (closure design does not permit the use of Cryocolourcode)

Catalogue No	Alt. No	Bottom	Volume, working, mL	Writing area [Y/ N]	Length, mm	Length [max.], mm	Pack qty
⑯ CRY-960-255M	375299	U	1.0	Yes	30	30	500
⑰ CRY-960-257Y	340711	U	1.8	Yes	44	44	500

Accessory

Catalogue No	Alt. No	Description	Dimensions, mm	Inner pack qty	Pack qty
⑱ CRY-960-505P	343958	Cover tubing, Nunc CryoFlex™, polyethylene	500 x 14	15	300

U = Round V = Conical