

Cylinders Class B



Cylinders, borosilicate glass, spouted, hexagonal base, heavy duty, Class B, DIN 12680, BS 604, ISO 4788

PYREX®

01

White enamel graduations. Thicker walls and heavily beaded rim for increased mechanical strength.

Catalogue No	Alt. No	Capacity, mL	Graduations, mL	Pack qty
CYP-340-070H	3280-50	50	1.0	2
CYP-340-090B	3280-100	100	1.0	2
CYP-340-110V	3280-250	250	2.0	2

Catalogue No	Alt. No	Capacity, mL	Graduations, mL	Pack qty
CYP-340-130P	3280-500	500	5.0	2
CYP-340-150J	3280-1000	1,000	10.0	2
CYP-340-170D	3280-2000	2,000	20.0	1



Cylinders, borosilicate glass, spouted, hexagonal base, heavy duty, squat form, Class B, DIN 12680, BS 604, ISO 4788

PYREX®

01

White enamel graduations. Thicker walls and heavily beaded rim for increased mechanical strength.

Catalogue No	Alt. No	Capacity, mL	Graduations, mL	Pack qty
CYP-360-050W	3285-25	25	1.0	2
CYP-360-070Q	3285-50	50	2.0	2
CYP-360-090K	3285-100	100	2.0	2
CYP-360-110H	3285-250	250	10.0	2

Catalogue No	Alt. No	Capacity, mL	Graduations, mL	Pack qty
CYP-360-130B	3285-500	500	10.0	2
CYP-360-150S	3285-1000	1,000	20.0	2



Cylinders, borosilicate glass, spouted, detachable base, Class B, **PYREX®**, DIN 12680, BS 604, ISO 4788

01

Glass cylinder with detachable plastic base. White graduations.

Catalogue No	Alt. No	Capacity, mL	Graduations, mL	Pack qty
CYP-330-030D	3205/04MP	10	0.2	2
CYP-330-050U	3205/06MP	25	0.5	2
CYP-330-070X	3205/08P	50	1.0	2
CYP-330-090Y	3205/10P	100	1.0	2
CYP-330-110F	3205/12P	250	2.0	2

Plastic bases only

Catalogue No	Alt. No	Capacity, mL	Graduations, mL	Pack qty
CYP-415-031A	4500/00	10	-	10
CYP-415-051R	4500/01	25	-	10
CYP-415-071F	4500/02M	50	-	10
CYP-415-091C	4500/04M	100	-	10
CYP-415-111L	4500/06M	250	-	10

Glass cylinders only

Catalogue No	Alt. No	Capacity, mL	Graduations, mL	Pack qty
CYP-371-030P	3205/05MP	10	0.2	2
CYP-371-050J	3205/07MP	25	0.5	2
CYP-371-070D	3205/09P	50	1.0	2
CYP-371-090U	3205/11P	100	1.0	2
CYP-371-110R	3205/13P	250	2.0	2