

Cylinders, borosilicate glass, polyethylene stopper and foot, Class B, DIN 12680, BS 604, ISO 4788

PYREX®

01

With robust plastic base and stopper. White graduations.

Catalogue No	Alt. No	Capacity, mL	Graduations, mL	Pack qty
CYP-410-050D	3215/06MP	25	0.5	2
CYP-410-070U	3215/08P	50	1.0	2
CYP-410-090X	3215/10P	100	1.0	2
CYP-410-110L	3215/12P	250	2.0	2

Polyethylene bases

Catalogue No	Alt. No	Capacity, mL	Pack qty
CYP-415-051R	4500/01	25	10
CYP-415-071F	4500/02M	50	10
CYP-415-091C	4500/04M	100	10
CYP-415-111L	4500/06M	250	10



Cylinders, borosilicate glass, hollow glass stopper, hexagonal base, Class B, DIN 12680, BS 604, ISO 4788

PYREX®

01

With glass stopper for chemical inertness. White graduations.

Catalogue No	Alt. No	Capacity, mL	Graduations, mL	Pack qty
CYQ-430-030Y	3211/04MP	10	0.2	2
CYQ-430-050C	3211/06MP	25	0.5	2
CYQ-430-070T	3211/08MP	50	1.0	2
CYQ-430-090N	3211/10MP	100	1.0	2
CYQ-430-130E	3211/14MP	500	5.0	2
CYQ-430-150V	3211/16MP	1,000	10.0	2
CYQ-430-170P	3211/18MP	2,000	20.0	1



Cylinders, Duran® glass, spouted, hexagonal base, Class B, DIN 12680, BS 604, ISO 4788



01

White graduations.

Catalogue No	Alt. No	Capacity, mL	Graduations, mL	Pack qty
CYL-520-010R	213960707	5	0.1	10
CYL-520-030L	213960801	10	0.2	10
CYL-520-050F	213961403	25	0.5	10
CYL-520-070W	213961703	50	1.0	10
CYL-520-090Q	213962408	100	1.0	10
CYL-520-110N	213963601	250	2.0	10
CYL-520-130H	213964409	500	5.0	10
CYL-520-150B	213965405	1,000	10.0	1
CYL-520-170S	213966307	2,000	20.0	1

