

Extraction thimbles, glass fibre

69

Whatman®

NEW

- Smoke stack gas monitoring
- Soxhlet extraction
- Analysing pesticide residues
- Determining oil/fat content of foods (e.g. French fries)
- Analysis of oil and grease in solid wastes
- Available in a range of sizes and wall thicknesses to suit your application
- Designed to fit most commercially available Soxhlet extractors

High-Purity Glass Microfibre Thimbles - The thimbles are completely free of binders or additives and can be used at temperatures up to 500°C or when using solvents incompatible with cellulose thimbles.

These thimbles are also used in pollution monitoring techniques (0.8µm nominal particle retention).

Typical thickness 1.7 mm.

Standard glass microfibre extraction thimbles

Grade 603g (Borosilicate glass fibre with inorganic binder)

Catalogue No	Alt. No	Size, mm	Thickness, mm	Pack qty
EXP-801-005K	10371005	16 x 50	1.0	25
EXP-801-010R	10371007	19 x 90	1.0	25
EXP-801-015H	10371011	22 x 80	1.5	25
EXP-801-020X	10371019	25 x 100	1.5	25
EXP-801-025E	10371023	26 x 100	1.5	25
EXP-801-030L	10371025	28 x 60	1.5	25
EXP-801-035B	10371029	25 x 98	1.5	25
EXP-801-040Y	10371036	30 x 100	1.5	25
EXP-801-045V	10371042	33 x 94	1.5	25
EXP-801-050F	10371043	33 x 100	1.5	25
EXP-801-055S	10371045	33 x 118	1.5	25
EXP-801-060C	10371055	35 x 150	1.5	25
EXP-801-065P	10371075	44 x 230	1.5	25
EXP-801-070W	10371103	10 x 38	1.5	25
EXP-801-075M	10371114	26.8 x 68	1.5	25
EXP-801-080T	10371122	26 x 63	1.5	25

High-Purity glass or quartz microfibre

Glass microfibre thimbles - Grade HP-GF

Made from high purity quartz microfibre, this thimble is able to withstand high temperatures (up to 1,000°C). Suitable for both solvent extraction and air sampling applications.

Catalogue No	Alt. No	Size, mm	Pack qty
EXP-801-205C	2814-199	19 x 90	25
EXP-801-210J	2812-259	Tapered, 25 x 90	10
EXP-801-215W	2814-300	30 x 100	25
EXP-801-220G	2814-432	43 x 123	25
EXP-801-225T	2814-533	53 x 135	25

Quartz microfibre (silica) thimbles

Catalogue No	Alt. No	Size, mm	Pack qty
EXP-801-230D	2814-259	Tapered, 25 x 90	25

