



# Personal Protective Equipment

## Hand Protection Disposable



Gloves, nitrile, disposable, blue, non-sterile, Semperguard, PPE category III



S3

Manufactured from high quality soft blue nitrile film, available in lightly powdered and powder free styles, both have a textured outer surface for better grip especially in wet conditions. Nitrile offers higher elasticity than latex, has higher puncture resistance and offers outstanding resistance against oils and fats. Ambidextrous with a rolled cuff, 0.25mm thickness at the palm and an AQL of 1.5. Longer cuff in sizes large and extra large.

Chemical resistance data for this glove is available at [www.sempermed.com](http://www.sempermed.com). Select the industrial gloves section from European Branch and then click on 'English', then 'Products', and then from the list headed 'Disposable Protective Gloves' select the appropriate nitrile glove and then click on 'Show List' for chemicals. Alternatively you can access the 'Products' section of the Sempermed site for this information via the Fisher Scientific website.

### Lightly powdered

Catalogue No	Size	Length, mm	Pack qty
SAR-255-030T	Small	240	100
SAR-255-050N	Medium	240	100
SAR-255-070H	Large	245	100

### Powder free

Catalogue No	Size	Length, mm	Pack qty
SAR-265-030M	Small	240	100
SAR-265-050G	Medium	240	100
SAR-265-070A	Large	245	100
SAR-265-090R	Extra large	255	90



Gloves, nitrile, disposable, powder free, Purple Nitrile™ and Purple Nitrile XTRA™, PPE Category III



S3

Thickness 0.11mm (palm) 0.15mm (textured fingerip). AQL 0.65.

### Purple Nitrile™, length 242mm

Catalogue No	Size	Pack qty
SAR-690-011Q	Extra small	100
SAR-690-031K	Small	100
SAR-690-051E	Medium	100
SAR-690-071V	Large	100
SAR-690-091P	Extra large	90

### Purple Nitrile XTRA™, length 305mm

Catalogue No	Size	Pack qty
SAR-692-011E	Extra small	50
SAR-692-031V	Small	50
SAR-692-051P	Medium	50
SAR-692-071J	Large	50
SAR-692-091D	Extra large	50