



Gloves, nitrile, disposable, non-sterile Touch 'n' Tuff™, PPE Category III

S3

- Provides three times the puncture resistance of natural latex and vinyl gloves
- No additives to contaminate product or irritate skin
- Rolled cuff, 0.1mm thickness at the palm
- Length 240mm

Powder free, green

Catalogue No	Size	Pack qty
SAR-480-100R	Small	100
SAR-480-200N	Medium	100
SAR-480-300J	Large	100
SAR-480-400F	Extra large	100


Powdered, green

Catalogue No	Size	Pack qty
SAR-480-010S	Small	100
SAR-480-020P	Medium	100
SAR-480-030M	Large	100
SAR-480-040J	Extra large	100

Powder free, blue with textured fingertips

Catalogue No	Size	Pack qty
SAR-605-100E	Small	100
SAR-605-200A	Medium	100
SAR-605-300T	Large	100
SAR-605-400P	Extra large	100

Ansell



Gloves, nitrile, disposable, non-sterile, powder free, Touch 'n' Tuff, green, PPE Category III

S3

Provides three times the puncture resistance of natural latex and vinyl gloves

- Length 300mm
- Exceptionally soft feel
- Proven splash resistance against hazardous chemicals
- Rolled cuff, 0.1mm thickness at the palm
- EN 374

Catalogue No	Size	Pack qty
SAR-491-020C	Small	100
SAR-491-030W	Medium	100
SAR-491-040T	Large	100
SAR-491-050Q	Extra large	100

Ansell



Gloves, nitrile, disposable, non-sterile, TNT, blue, PPE Category III

S3

NEW

- Provides three times the puncture resistance of natural latex and vinyl gloves
- Exceptionally soft feel
- No additives to contaminate product or irritate skin
- Rolled cuff, 0.1mm thickness at the palm
- Length 300mm

Catalogue No	Size	Pack qty
SAR-600-010M	Small	100
SAR-600-020J	Medium	100
SAR-600-030G	Large	100
SAR-600-040D	Extra large	100

Ansell

