



Disposable flasks, Erlenmeyer, PETG, sterile, Thermo Scientific Nalgene®

Thermo
SCIENTIFIC

54



Optically clear PETG plastic. Ideal for shaker and suspension cell culture, media preparation and storage. Gamma sterilised, moulded-in graduations. Supplied with HDPE screw closures.

Plain bottom, white closure

Catalogue No	Alt. No	Capacity, mL	Graduations, mL	Pack qty
TKV-137-005E	4112-0125	125	25	24
TKV-137-010L	4112-0250	250	50	12
TKV-137-020Y	4112-0500	500	50	12
TKV-137-025V	4112-1000	1,000	50	6
TKV-137-030F	4112-2000	2,000	500	4
TKV-137-030W	4112-2800	2,800	250	4

Baffled bottom, white closure

Catalogue No	Alt. No	Capacity, mL	Graduations, mL	Pack qty
TKV-137-035S	4112-0125	125	25	24
TKV-137-040C	4113-0250	250	50	12
TKV-137-050W	4113-0500	500	50	12
TKV-137-055M	4113-1000	1,000	50	6
TKV-137-060T	4113-2000	2,000	500	4
TKV-137-065J	4113-2800	2,800	250	4

Plain bottom, blue vented closure

Catalogue No	Alt. No	Capacity, mL	Graduations, mL	Pack qty
TKV-137-200G	4115-0125	125	25	24
TKV-137-210D	4115-0250	250	50	12
TKV-137-220A	4115-0500	500	50	12
TKV-137-230U	4115-1000	1,000	50	6
TKV-137-240R	4115-2000	2,000	500	4
TKV-137-245H	4115-2800	2,800	250	4

Baffled bottom, blue vented closure

Catalogue No	Alt. No	Capacity, mL	Graduations, mL	Pack qty
TKV-137-300C	4116-0125	125	25	24
TKV-137-310W	4116-0250	250	50	12
TKV-137-320T	4116-0500	500	50	12
TKV-137-330Q	4116-1000	1,000	50	6
TKV-137-340N	4116-2000	2,000	200	4
TKV-137-350K	4116-2800	2,800	250	4

Vented closure, 0.2mm hydrophobic PTFE membrane

Catalogue No	Alt. No	Description	Pack qty
TKV-137-100K	4114-0038	Fits 250mL flasks	12
TKV-137-110H	4114-0045	Fits 500, 1,000 and 2,000mL flasks	12



Flasks, shaker, borosilicate glass, three baffles

03

WHEATON
Science International



Designed for use on either rotary or reciprocating shakers for cell or bacterial cultures.

Special design with three baffles maximises aeration.

Provided with break-resistant lip and marking spot on side of flask. Break-resistant lip accommodates cotton plug or stopper. Wheaton-33® low extractable borosilicate glass (meets ASTM Type I Class A and USP Type I standards).

Catalogue No	Alt. No	Volume, mL
FHB-450-030T	353255	50
FHB-450-050N	353257	125
FHB-450-070H	353259	250
FHB-450-090B	353262	500
FHB-450-110V	353264	1,000
FHB-450-130P	353266	2,000