



Flasks, Pyrex® glass, Kjeldahl

01

PYREX®

Long neck. Suitable for determination of nitrogen.

Catalogue No	Alt. No	Capacity, mL	Neck I.D., mm	Height, mm
FHD-430-030N	1200/04M	100	17	240

Catalogue No	Alt. No	Capacity, mL	Neck I.D., mm	Height, mm
FHD-430-050H	1200/10D	300	24	300

Catalogue No	Alt. No	Capacity, mL	Neck I.D., mm	Height, mm
FHD-430-070B	1200/12D	500	30	325
FHD-430-090S	1200/16M	800	30	350



Flasks, borosilicate glass, iodine, BS 2735

01



Catalogue No	Alt. No	Capacity, mL	Pack qty
QFY-952-E	FIB 250/3	250	2
QFY-972-N	FIB 500/3	500	2

Iodine stopper

Catalogue No	Alt. No	Pack qty
QTE-722-S	SIB24	2



Flasks, polypropylene, conical, filter, Thermo Scientific Nalgene®

55

Thermo
SCIENTIFIC


Moulded in one piece, will withstand full vacuum. Tubulure at 45° angle prevents tipping. Adapter accepts 6.4mm to 9.5mm I.D. tubing. Autoclavable.

Catalogue No	Alt. No	Capacity, mL	Neck I.D., mm	Pack qty
FHB-580-010T	DS4101-0500	500	28	1
FHB-580-020Q	DS4101-1000	1,000	31	1



Flasks, PMP, conical, Thermo Scientific Nalgene®

55

Thermo
SCIENTIFIC


Erlenmeyer, transparent, with good chemical resistance. Graduated - with polypropylene screw closure. Autoclavable.

Catalogue No	Alt. No	Capacity, mL	Pack qty
FHB-690-040S	4109-0125	125	6
FHB-690-050P	4109-0250	250	4
FHB-690-060J	4109-0500	500	4
FHB-690-110A	4109-1000	1,000	2