



Personal Protective Equipment

Hand Protection Liners



Glove liners, cleanroom

03



100% continuous filament knitted polyester, designed to be worn under latex, PVC or other outer gloves for maximum comfort and dexterity. One size. May be cleanroom laundered.

Catalogue No	Size	Type	Pack type	Pack qty
SAR-640-010H	Standard	Half-fingered	Pair	10
SAR-640-020E	Large	Half-fingered	Pair	10
SAR-640-040V	Standard	Full-fingered	Pair	10
SAR-640-050S	Large	Full-fingered	Pair	10



Glove liners, thermal insulating Therm-A-Knit™

S4



Light, stretchable thermal protection glove for all cold conditions, made from DuPont™ ThermoStat™ fabric. Hollow core fibres for greater insulation and warmth. Rapid perspiration elimination through fibres maintains body heat. Can be worn on their own or as liners, particularly in the laboratory beneath nitrile examination gloves where work is being undertaken in cold and damp conditions. Ambidextrous. No seams to rub or chafe. Washable at maximum 40°C. Colour blue.

Catalogue No	Size	Length, mm	Pack type	Pack qty
SAR-572-030W	7	240	Pair	12
SAR-572-070K	9	260	Pair	12



Glove liners, Insulator™

S4



- Superior thermal insulation for warmth
- Manufactured using Thermax™ fibre from DuPont™
- Soft and non-chaffing with a fast wicking action for greater dryness
- Machine washable/dryable
- Use as a liner for industrial rubber gloves or on its own
- Ideal for food and fishing industries, construction, food retailing and general outdoor maintenance work in cold weather
- For minimal risks only

Catalogue No	Pack type	Pack qty
SAR-565-010R	Pair	12



Gloves, polyester, Monysoft

S4

For product protection, assembly operations or inspection tasks.

Lint free glove with excellent comfort and dexterity. Seamless construction for comfort.

Catalogue No	Size	Pack qty
SAR-863-010F	6 to 7	1
SAR-863-020C	7 to 8	1
SAR-863-030W	8 to 9	1
SAR-863-040T	9 to 10	1