

Flasks, borosilicate glass, volumetric, Class A, ISO 1042, DIN 12664

Fisherbrand®

41

With one mark graduation and interchangeable polyethylene stopper. The details on each flask include confirmation of compliance to International Standards as well as an individual re-order code.

Clear, white markings

Catalogue No	Volume, mL	Tolerance, mL	Stopper diameter, mm	Pack qty
FB50141	5	0.025	10/19	5
FB50143	10	0.025	10/19	5
FB50145	20	0.04	10/19	5
FB50147	25	0.04	10/19	5
FB50149	50	0.06	12/21	5
FB50151	100	0.10	12/21	5
FB50153	200	0.15	14/23	2
FB50155	250	0.15	14/23	2
FB50157	500	0.25	19/24	2
FB50159	1,000	0.40	24/29	2
FB50161	2,000	0.60	29/32	1

Amber, white markings

Catalogue No	Volume, mL	Tolerance, mL	Stopper diameter, mm	Pack qty
FB55531	5	0.025	10/19	5
FB55533	10	0.025	10/19	5
FB55535	20	0.04	10/19	5
FB55537	25	0.04	10/19	5
FB55539	50	0.06	12/21	5
FB55541	100	0.10	12/21	5
FB55543	200	0.15	14/23	2
FB55545	250	0.15	14/23	2
FB55547	500	0.25	19/24	2
FB55549	1,000	0.40	24/29	2



Flasks, borosilicate glass, volumetric, trapezoidal, Class A

volac

01

Mini volumetric flasks with tight accuracy tolerances, and shaped to maximise stability while allowing access to contents. Individually calibrated and tested.

Clear

Catalogue No	Alt. No	Capacity, mL	Tolerance, mL	Pack qty
FHM-440-010X	RB257/A/5	1	0.010	5
FHM-440-030Y	RB257/B/5	2	0.015	5
FHM-440-050C	RB257/C/5	3	0.025	5
FHM-440-070T	RB257/D/5	4	0.025	5
FHM-440-090N	RB257/E/5	5	0.025	5

Amber

Catalogue No	Alt. No	Capacity, mL	Tolerance, mL	Pack qty
FHM-450-010H	RB257/AM/A/5	1	0.010	5
FHM-450-030B	RB257/AM/B/5	2	0.015	5
FHM-450-050S	RB257/AM/C/5	3	0.025	5
FHM-450-070M	RB257/AM/D/5	4	0.025	5
FHM-450-090G	RB257/AM/E/5	5	0.025	5

