

Funnels glass conical



Funnels, soda lime glass, conical

41

Fisherbrand®

Catalogue No	Top diameter, mm	Stem length, mm	Pack qty
FB56630	35	40	1
FB56631	45	45	1
FB56632	50	50	1
FB56633	55	55	1
FB56634	60	60	1
FB56635	70	70	1
FB56636	100	100	1



Funnels, borosilicate glass, conical

41

Fisherbrand®

Catalogue No	Top diameter, mm	Pack qty
FB56640	55	10
FB56641	75	10
FB56642	100	10
FB56643	150	1
FB56644	200	1



Funnels, Pyrex® glass, conical

01

PYREX®

With short stem.

Catalogue No	Alt. No	Top diameter, mm	Stem length, mm	Stem O.D., mm
FPH-360-005J	2140/01M	35	35	8
FPH-360-010Q	2140/02M	45	45	8
FPH-360-020N	2140/03M	65	55	8
FPH-360-030K	2140/06M	75	75	8
FPH-360-050E	2140/10M	100	100	10
FPH-360-060B	2140/11M	125	125	12
FPH-360-070V	2140/12M	150	150	16
FPH-360-090P	2140/14M	200	175	24



Funnels, Duran® glass, conical, DIN 12445, ISO 4798

01



With short stem.

Catalogue No	Alt. No	Top diameter, mm	Stem length, mm	Stem O.D., mm	Filter paper diameter, mm	Pack qty
FPH-310-010F	213512308	35	35	6	45 to 55	10
FPH-310-030W	213512805	45	45	6	55 to 70	10
FPH-310-050Q	213513304	55	55	8	70 to 90	10
FPH-310-070K	213513801	70	70	8	110 to 125	10
FPH-310-090E	213514103	80	80	10	125 to 150	10
FPH-310-110B	213514609	100	100	10	150 to 185	10
FPH-310-130S	213515708	150	150	16	240 to 270	10
FPH-310-150M	213516104	200	175	26	320 to 385	1
FPH-310-170G	213516901	300	175	30	500	1