

Funnels plastic conical



Funnels, polypropylene, conical, powder

02

Kartell



60° angled sides with external ribs avoid pressure lock. Large bore avoids powder bridging during transfer. Autoclavable.

Catalogue No	Top diameter, mm	Stem O.D., mm	I.D., mm	Stem length, mm	Pack qty
FPH-750-007D	60	15	11	18	5
FPH-750-011M	80	15	11	22	5
FPH-750-021J	100	25	22	23	5
FPH-750-030Y	120	30	26	27	1
FPH-750-036T	150	35	33	41	1
FPH-750-040F	180	43	40	47	1



Funnels, PMP, conical, filter

02



Unsurpassed for non-vacuum, high-speed filtration, these clear funnels have internal helicoid ribs to increase the filtration speed of filter paper. Vertical outer ribs prevent airlock. Unaffected by corrosive liquids and heat resistant to 177°C (350°F).

Catalogue No	Alt. No	Top diameter, mm	I.D., mm	Stem O.D., mm	Stem length, mm	Filter paper diameter, mm	Pack qty
FPH-430-030A	F14640-0000	51	3	7.5	150	75/90	6
FPH-430-050R	F14642-0000	70	3	7.5	150	110/125	6
FPH-430-070L	F14644-0000	70	7	10.0	70	110/125	6
FPH-430-090F	F14646-0000	100	7	10.0	108	150/185	4
FPH-430-110C	F14648-0000	140	10	16.0	132	240/270	3
FPH-430-130T	F14650-0000	196	20	26.0	160	320/385	2



Funnels, LDPE, conical, Thermo Scientific Nalgene®

55



Large funnels with tabs.

Catalogue No	Alt. No	Capacity, mL	Top diameter, mm	Stem O.D., mm	Material	Pack qty
FPH-465-010B	4260-0020	29	52	10	LDPE	6
FPH-465-020V	4260-0030	104	76	11	LDPE	6
FPH-465-030S	4260-0040	254	101	12	LDPE	2
FPH-465-040P	4260-0060	777	152	15	LDPE	1
FPH-465-050M	4260-0100	3,684	248	23	LDPE	1
FPH-465-060J	4260-0120	6,976	343	33	HDPE	1



Funnels, PTFE

41



Chemically inert with super smooth, non-stick internal finish.

Catalogue No	Top diameter, mm	Stem O.D., mm	Height, mm
FB58870	30	8	50
FB58871	50	10	85
FB58872	70	12	110
FB58873	100	14	170
FB58874	150	16	220