

Dispersers T 18 basic and T 25 digital ULTRA-TURRAX® 04



2 yr

Dispersing elements, stainless steel, for T 25 digital

Catalogue No	Alt. No	Description
TWK-610-010M	S 25N-8G	Working volume 1 to 50mL, 108mm length x 8mm dia., rotor dia. 6.1mm, min./max. immersion depth 27/85mm
TWK-610-035M	S 25N-10G	Working volume 1 to 100mL, 105mm length x 10mm dia., rotor dia. 7.5mm, min./max. immersion depth 22/85mm
TWK-610-045J	S 25N-18G	Working volume 10 to 1,500mL, 194mm length x 18mm dia., rotor dia. 12.7mm, min./max. immersion depth 40/165mm
TWK-610-050A	S 25N-25G	Working volume 50 to 2,000mL, 194mm length x 25mm dia., rotor dia. 17mm, min./max. immersion depth 40/165mm
TWK-610-070R	S 25N-25F	Working volume 100 to 2,000mL, 194mm length x 25mm dia., rotor dia. 18mm, min./max. immersion depth 40/165mm

Dispersing elements, plastic, FDA approved, autoclavable, for T 18 basic
(Further dispersing tools are available on request)

Catalogue No	Alt. No	Description	Pack qty
TWK-625-500U	S18D-10G-KS	Working volume 10 to 100mL, 150mm length x 10mm dia., rotor dia. 6.75mm, min./max. immersion depth 15/85mm	5
TWK-625-505K	S18D-10G-KS	As above	10
TWK-625-510R	S18D-14G-KS	Working volume 10 to 500mL, 150mm length x 14mm dia., rotor dia. 9.5mm, min./max. immersion depth 15/85mm	5
TWK-625-515H	S18D-14G-KS	As above	10

Dispersing elements, plastic, FDA approved, autoclavable, for T 25 digital
(Further dispersing tools are available on request)

Catalogue No	Alt. No	Description	Pack qty
TWK-625-520X	S25D-10G-KS	Working volume 10 to 100mL, 150mm length x 10mm dia., rotor dia. 6.75mm, min./max. immersion depth 15/85mm	5
TWK-625-525E	S25D-10G-KS	As above	10
TWK-625-530L	S25D-14G-KS	Working volume 10 to 500mL, 150mm length x 14mm dia., rotor dia. 9.5mm, min./max. immersion depth 15/85mm	5
TWK-625-535B	S25D-14G-KS	As above	10

Accessories and spares

Catalogue No	Alt. No	Description
SWR-535-515T	3160000	Plate stand, 200mm x 316mm x 560mm (w x d x h)
TWK-614-505A	2657700	Bosshead
TWK-614-510H	3008600	Strap clamp
TWK-614-508R	DZM control. O	Revolution counter



*Polymerase Chain Reaction (PCR) is a process covered by patents owned by Hoffman-La Roche

Are you online?

- See the price you pay
- Live stock availability
- Order tracking
- Fully secure online ordering

Visit our website

