


Inoculating

Inoculating loops, polystyrene, sterile

S



High impact, polystyrene, disposable inoculating loops. 200mm long to reduce the risk of contamination.


Catalogue No	Alt. No	Capacity, μL	Inner pack qty	Pack qty
LPS-131-011B	PL2010/1000	10	20	1,000
LPS-131-021V	PL205/1000	5	20	1,000
LPS-131-031S	PL201/1000	1	20	1,000

Inoculating loops, with sphere, sterile, Quadloops®

56



S



Sterile disposable inoculating device for either sampling fixed volumes or for dilution streaking to obtain isolated colonies.

Unique sphere/handle form provides up to 4 sterile streaking surfaces simply by turning the handle 90° after each streaking, thereby exposing a new sterile streaking surface.

Catalogue No	Alt. No	Capacity, μL	Colour	Pack qty
DIS-170-015G	QL1	1	Green	1,000
DIS-170-025D	QL10	10	Blue	1,000

03

Inoculating needles, polystyrene

03



S



Clear polystyrene with chisel tip to facilitate stab cultures and recovery of adherent colonies.

Catalogue No	Description	Inner pack qty	Pack qty
DIS-107-010A	Clear polystyrene with chisel tip	20	1,000
DIS-107-030R	Clear polystyrene with chisel tip individually wrapped	1	1,000

Inoculating loops, sterile

03



S



Manufactured to the highest quality using cadmium free pigments.

10 μL green and 1 μL clear loops are manufactured in hard crystal polystyrene for general plating. 10 μL white, 5 μL red and 1 μL yellow are manufactured out of soft HI polystyrene, making them ideal for clinical use as there is no danger of the head breaking off in use.

Catalogue No	Colour	Capacity, μL	Inner pack qty	Pack qty
DIS-105-010M	White	10	20	1,000
DIS-105-020J	Green	10	20	1,000
DIS-105-030G	Red	5	20	1,000
DIS-105-040D	Clear	1	20	1,000
DIS-105-050A	Yellow	1	20	1,000

Individually wrapped

Catalogue No	Colour	Capacity, μL	Pack qty
DIS-105-015C	White	10	1,000
DIS-105-025W	Green	10	1,000
DIS-105-035T	Red	5	1,000
DIS-105-045Q	Clear	1	1,000
DIS-105-055N	Yellow	1	1,000