

## Label printer, thermal transfer, LAB PAL™

50



- Easy to use, hand held thermal transfer printer
- Includes time/date stamp, banner mode, serialisation font sizes
- Brady's labels resist most chemicals - xylene, alcohol, ethanol, DMSO

The labels will stick and stay stuck even after exposure to harsh environments including liquid nitrogen, freezers and hot water baths as well as chemicals.

The drop-in cartridge makes loading quick and the exposed edge on the labels makes peeling easy, even with gloves. It is lightweight and runs on 6 x AA batteries. To be applied at room temperature.

### Technical Specification

Font size, point	7, 10, 14, 20 and 28
Operating temperature, °C	4 to 49
Dimensions [l x w x h], mm	228.6 x 114.3 x 63.5

### Applications key:

Environmental conditions	Product identification
F Freezer (-70°C)	Low temperature freezer
N Liquid nitrogen (-196°C)	S Slide Identification
A Autoclave	P Plate Identification
H Hot water baths (100°C)	G General Identification
L Low-profile, (extremely thin)	

Label material: continuous tape width x 5.8m (approx)

Catalogue No	
LAB-932-010P	

### Accessories

Catalogue No	Description
LAB-932-500V	LABPAL adapter IDPAL-AC
BMT-680-050T	Battery, alkaline, 1.5V, AA

Label material: polyester. Applications: F N A H L / V P. Must overlap onto itself. Ideal for straws

Catalogue No	Alt. No	Size, mm
LAB-932-100X	PAL-375-461	9.5 x 6m
LAB-932-106C	PAL-500-461	12.7 x 6m
LAB-932-112H	PAL-750-461	19.1 x 6m

Label material: nylon cloth. Applications: F N A V / P G. Ideal for tubes, vials and bottles

Catalogue No	Alt. No	Size, mm
LAB-932-101M	PAL-375-499	9.5 x 5.8mm
LAB-932-107A	PAL-500-499	12.7 x 5.8mm
LAB-932-113F	PAL-750-499	19.1 x 5.8mm

Label material: polyester. Applications: F A S / P G. Ideal for slides, DMSO and xylene resistant

Catalogue No	Alt. No	Size, mm
LAB-932-102K	PAL-375-488	9.5 x 6m
LAB-932-108V	PAL-500-488	12.7 x 6m
LAB-932-114D	PAI-750-488	19.1 x 6m

Label material: polyester. Applications: F A / P G. Ideal for plates, bottles and flat surfaces

Catalogue No	Alt. No	Description	Size, mm
LAB-932-103Y	PAL-375-423	B423	9.5 x 6mm
LAB-932-104G	PAL-375-430	B430	9.5 x 6mm
LAB-932-109T	PAL-500-423	B423	12.7 x 6mm
LAB-932-110L	PAL-500-423	B430	12.7 x 6mm
LAB-932-115B	PAL-750-423	B423	19.1 x 6mm
LAB-932-116W	PAL-750-430	B430	19.1 x 6mm

Label material: vinyl. Applications: F / P G. Ideal for labelling racks

Catalogue No	Alt. No	Colour	Size, mm
LAB-932-105E	PAL-375-439	White	9.5 x 6m
LAB-932-111J	PAL-500-439	White	12.7 x 6m
LAB-932-117U	PAL-750-439	White	19.1 x 6m
LAB-932-120Y	PAI-375-439-BL	Blue	9.5 x 6m
LAB-932-121G	PAL-375-439-GR	Green	9.5 x 6m
LAB-932-122E	PAL-375-439-OR	Orange	9.5 x 6m
LAB-932-123C	PAL-375-439-RD	Red	9.5 x 6m
LAB-932-124A	PAL-375-439-YL	Yellow	9.5 x 6m

