


Microplates, 96 well, polystyrene,
Thermo Scientific Nunc®, MicroWell™, Nunclon™ Δ

Thermo
SCIENTIFIC

76

S

- Plate edges and lids designed to reduce evaporation
- Non-toxic
- Raised well rims reduce cross contamination risk
- Lids shaped to facilitate stacking and handling
- Excellent optical quality
- Certified
- Nunclon™ Δ surface for cell based assays



Catalogue No	Alt. No	Well shape	Lid	Volume, working, µL	Culture area, cm²	Inner pack qty	Pack qty
TKT-180-070U	167008	Flat	Yes	200	0.33/well	1	50
TKT-180-090X	168055	Flat	No	200	0.33/well	1	50
TKT-182-030U	161093	Flat	Yes	200	0.33/well	10	160
TKT-182-050X	156545	Flat	No	200	0.33/well	10	180
TKT-180-010P	143761	Round	No	250	0.36/well	1	50
TKT-180-050D	163320	Round	Yes	250	0.36/well	1	50
DIS-984-112F	168136	Round	Yes	250	0.36/well	10	160

Accessories - Lids for microplates

Catalogue No	Alt. No	Description	Radiation sterilised [Y/ N]	Inner pack qty	Pack qty
TKT-180-501T	263339	Lid	N	5	100
TKT-180-503P	264122	Lid	Y	1	50
TKT-180-502R	264623	Lid	Y	20	60

Microplates, 96 well, polystyrene, optical bottom,
Thermo Scientific Nunc®, MicroWell™, Nunclon™ Δ

Thermo
SCIENTIFIC

76

S CE IVD

- Glass bottom plates are with No. 1.5 German cover glass for minimum light scatter and low autofluorescence
- Polymer bottom plates combine optical clarity of virgin crystalline polystyrene with optimal surface for wide range of HTS applications
- Flat bottom well geometry for plate reader access, working volume 200µL, culture area 0.32cm²/well
- Polystyrene upper structure available in black, white or clear
- TC-treated



Note: Glass thickness No. 1 = 0.13 to 0.16mm. No. 1.5 = 0.16 to 0.19mm.

No.1 cover glass plates

Catalogue No	Alt. No	Upper colour structure	Lid	Inner pack qty	Pack qty
TKT-195-020V	164588	Black	Yes	6	30
TKT-195-030S	164590	White	Yes	6	30

Polymer bottom plates

Catalogue No	Alt. No	Upper colour structure	Lid	Inner pack qty	Pack qty
TKT-195-040P	165306	White	Yes	10	30
TKT-195-060J	165305	Black	No	10	30

Lids, clear

Catalogue No	Alt. No	Sterile, Yes/No	Inner pack qty	Pack qty
TKT-180-501T	263339	No	5	100
TKT-180-503P	264122	Yes	1	50
MPA-740-504F	250002	No	20	60
MPA-740-506B	250003	No	20	60
MPA-740-510K	250005	Yes	20	180