



## Indicator paper, sticks, pH, non-bleeding

41 **Fisherbrand®**

85mm x 6mm sticks with different indicator papers sealed onto each strip.

The pH is determined by comparison with a colour chart supplied. The sticks are sufficiently long to avoid contact between the fingers and the test sample.

Catalogue No	Range, pH	Pack qty
FB33003	0 to 14	100
FB33005	0 to 6	100
FB33007	1.7 to 3.8	100
FB33009	3.6 to 6.1	100
FB33011	4.5 to 10	100
FB33013	6.0 to 7.7	100
FB33015	7.0 to 14	100



## Indicator paper, reels

41 **Fisherbrand®**

Dispenser reels each, 5m x 7mm (l x w), except FB33057 which is 5m x 10mm (l x w).

Catalogue No	Type	Range, pH	Colour change	Pack qty
FB33033	Litmus red	5 to 8	Red to blue	4
FB33035	Litmus blue	5 to 8	Blue to red	4
FB33037	Congo red	3 to 5	Red to blue	4
FB33039	Phenolphthalein	8.5 to 10	Colourless to red	4
FB33041	pH	1 to 14	-	4
FB33043	pH	1 to 11	-	4
FB33045	pH	0.5 to 5.5	-	4
FB33047	pH	4.0 to 7.0	-	4
FB33049	pH	6.4 to 8.0	-	4
FB33051	pH	8 to 10	-	4
FB33053	pH	9 to 13	-	4
FB33055	pH	12 to 14	-	4
FB33057	pH	1 to 11	-	4



## Indicator paper

69 **Whatman®**

### Strips

**Type CF:** 80mm x 6mm (l x w), each with four individual segments of different dye-impregnated paper mounted on plastic support.

**Type CS:** 100mm x 11mm (l x w). Every test strip has its own printed colour comparison chart of eight or more marked pH values, with a central indicator segment.

### Reels

Three colour bands giving unique colour combinations for each test. These are compared with the eleven colour coded pH values in the chart printed on the dispenser insert.

### Strips

Catalogue No	Alt. No	Range, pH	Type	Pack qty
PBK-380-010G	2613-991	0 to 14	CF	100
PBK-380-020D	2612-990	1 to 12	CS	200
PBK-380-030A	2626-990	1.8 to 3.8	CS	200
PBK-380-040U	2627-990	3.8 to 5.5	CS	200
PBK-380-050R	2628-990	5.2 to 6.8	CS	200
PBK-380-060X	2629-990	6.0 to 8.1	CS	200
PBK-380-070L	2630-990	8.0 to 9.7	CS	200

### Reels

Catalogue No	Alt. No	Range, pH	Dimensions	Pack qty
PBK-390-030Q	2611-628	1 to 11	5m x 10mm	1 reel

### Books

Catalogue No	Alt. No	Range, pH	Pack type	Pack qty
PBK-340-030F	2600-500	pH 1 to 11	10 books (10 strips per book)	20
PBK-340-050W	2600-601	Litmus Blue	10 books (20 leaves per book)	10
PBK-340-070Q	2600-602	Litmus Red	10 books (20 leaves per book)	10

Specialised test papers for H<sub>2</sub>S, chlorine and other oxidising reagents

Catalogue No	Alt. No	Type	Pack type	Pack qty
PBK-350-010E	2602-501A	Lead acetate	Dispenser - 7mm x 5m	1
PBK-350-020B	2602-500A	Potassium iodide	Dispenser - 7mm x 5m	1
PBK-350-025X	2651-500	Starch Iodine	Books of 20 strips	10